

# **ab175811 – 14,15-DHET ELISA Kit**

## Instructions for Use

A competitive immunoenzymatic assay for the quantitative measurement of 14,15-DHET in serum, plasma, urine, cell culture extracts and tissues.

This product is for research use only and is not intended for diagnostic use.

# Table of Contents

---

## INTRODUCTION

- |                  |   |
|------------------|---|
| 1. BACKGROUND    | 2 |
| 2. ASSAY SUMMARY | 3 |

## GENERAL INFORMATION

- |                                     |   |
|-------------------------------------|---|
| 3. PRECAUTIONS                      | 4 |
| 4. STORAGE AND STABILITY            | 4 |
| 5. MATERIALS SUPPLIED               | 4 |
| 6. MATERIALS REQUIRED, NOT SUPPLIED | 5 |
| 7. LIMITATIONS                      | 5 |
| 8. TECHNICAL HINTS                  | 6 |

## ASSAY PREPARATION

- |                                   |    |
|-----------------------------------|----|
| 9. REAGENT PREPARATION            | 7  |
| 10. STANDARD PREPARATION          | 8  |
| 11. SAMPLE COLLECTION AND STORAGE | 9  |
| 12. PLATE PREPARATION             | 15 |

## ASSAY PROCEDURE

- |                     |    |
|---------------------|----|
| 13. ASSAY PROCEDURE | 16 |
|---------------------|----|

## DATA ANALYSIS

- |                       |    |
|-----------------------|----|
| 14. CALCULATIONS      | 17 |
| 15. TYPICAL DATA      | 18 |
| 16. ASSAY SPECIFICITY | 19 |

## RESOURCES

- |                     |    |
|---------------------|----|
| 17. TROUBLESHOOTING | 20 |
| 18. NOTES           | 22 |

## 1. BACKGROUND


Abcam's 14,15-DHET competitive *in vitro* ELISA (Enzyme-Linked Immunosorbent Assay) kit is designed for accurate quantitative measurement in serum, plasma, urine, cell culture extracts and tissues following proper isolation and purification.

A 96-well plate has been precoated with 14,15-DHET antibody. Samples and the 14,15-DHET-HRP conjugate are added to the wells, where 14,15-DHET in the sample competes with the added 14,15-DHET-HRP for antibody binding. After incubation, the wells are washed to remove unbound material and TMB substrate is then added which is catalyzed by HRP to produce blue coloration. The reaction is terminated by addition of Stop Solution which stops the color development and produces a color change from blue to yellow. The intensity of signal is inversely proportional to the amount of 14,15-DHET in the sample and the intensity is measured at 450 nm.

The 14,15-DHET level exhibited strong positive correlation with hypertension, brain injury and stroke in rodents. 14,15-DHET is a representative metabolite of soluble epoxide hydrolase-mediated metabolism of EETs, which are generated by arachidonic acid epoxygenase activity of cytochromes P<sub>450</sub>. Human blood 14,15-DHET levels were measured using the 14,15-DHET ELISA kit. Increased 14,15-DHET levels of human cells as detected by the 14,15-DHET ELISA were indicative of the neoplastic and metastatic phenotype of carcinoma cells.


## 2. ASSAY SUMMARY

Capture Antibody




Prepare all reagents and samples as instructed.

Sample



Add standards and samples to each well used.

Labeled HRP-Conjugate



Add prepared HRP conjugate to each well and incubate at room temp.

Substrate **Colored Product**



Add TMB substrate to each well. Incubate at room temperature. Add Stop Solution to each well. Read immediately.

## 3. PRECAUTIONS

**Please read these instructions carefully prior to beginning the assay.**

All kit components have been formulated and quality control tested to function successfully as a kit. Modifications to the kit components or procedures may result in loss of performance.

## 4. STORAGE AND STABILITY

**Store kit at 2-8°C or -20°C immediately upon receipt.**

Refer to list of materials supplied for storage conditions of individual components. Observe the storage conditions for individual prepared components in section 9. Reagent Preparation.

## 5. MATERIALS SUPPLIED

Item	Amount	Storage Condition (After Preparation)
14,15-DHET ELISA plate	96 Wells	2-8°C
14,15-DHET Standard (1 mg/mL)	2 µL	2-8°C
1,000X 14,15-DHET HRP Conjugates	12 µL	2-8°C
10X Sample Dilution Buffer	25 mL	2-8°C
HRP Buffer	15 mL	2-8°C
10X Wash Buffer Solution	25 mL	2-8°C
TMB Substrate	24 mL	2-8°C

## 6. MATERIALS REQUIRED, NOT SUPPLIED

These materials are not included in the kit, but will be required to successfully utilize this assay:

- Microplate reader capable of measuring absorbance at 450 or 620 nm
- Incubator at 37°C
- Multi- and single-channel pipettes to deliver volumes between 10 and 1,000  $\mu\text{L}$
- Optional: Automatic plate washer for rinsing wells.
- Storage bottles
- Rotating mixer
- Deionised or (freshly) distilled water.
- Disposable tubes
- Timer

## 7. LIMITATIONS

- ELISA kit intended for research use only. Not for use in diagnostic procedures
- Use only clean pipette tips, dispensers, and lab ware
- Do not interchange screw caps of reagent vials to avoid cross-contamination
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate, without splashing, accurately to the bottom of wells

### 8. TECHNICAL HINTS

- Avoid foaming or bubbles when mixing or reconstituting components
- Avoid cross contamination of samples or reagents by changing tips between sample, standard and reagent additions
- Ensure plates are properly sealed or covered during incubation steps
- Complete removal of all solutions and buffers during wash steps is necessary for accurate measurement readings
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background
- **This kit is sold based on number of tests. A 'test' simply refers to a single assay well. The number of wells that contain sample, control or standard will vary by product. Review the protocol completely to confirm this kit meets your requirements. Please contact our Technical Support staff with any questions**

## 9. REAGENT PREPARATION

Equilibrate all reagents, samples and controls to room temperature (18-25°C) prior to use.

### 9.1 **1X Wash Buffer**

Mix the 10X Wash Buffer Solution with a stir bar, applying low, gentle heat until a clear colorless solution is obtained. Dilute the entire contents of the 10X Wash Buffer Solution (25 mL) with 225 mL of deionized water to yield a final volume of 250 mL of 1 X Wash Buffer. This can then be refrigerated for the entire life of the kit.

### 9.2 **1X HRP Conjugate**

Dilute 1 vial of the 14,15-DHET-HRP conjugate (12 µL) with 12 mL of HRP Buffer. One vial makes enough conjugate for one plate. The conjugate must be used the same day and should not be stored for later use.

### 9.3 **1X Sample Dilution Buffer**

Prepare 1X Sample Dilution Buffer by adding 25 mL of 10X Sample Dilution Buffer to 225 mL of dH<sub>2</sub>O. Mix gently and thoroughly.

## 10. STANDARD PREPARATION

Prepare serially diluted standards immediately prior to use. Always prepare a fresh set of positive controls for every use.

- 10.1 Label 5 microtubes as Standard # 2 - 6.
- 10.2 Add 900  $\mu$ L of the 1X Sample Dilution Buffer to the microtubes for Standards # 2 to 6.
- 10.3 Prepare a 1  $\mu$ g/mL **Standard #1** by first spinning down the enclosed 14,15-DHET standard vial (2  $\mu$ L, filled with inert gas) and then adding 1.998 mL of 1X Sample Dilution Buffer to obtain 2 mL of solution.
- 10.4 Prepare **Standard #2** by adding 100  $\mu$ L of the Standard #1 to the microtube labeled **Standard #2**. Mix thoroughly and gently.
- 10.5 Prepare **Standard #3** by adding 100  $\mu$ L of the Standard #2 to the microtube labeled **Standard #3**. Mix thoroughly and gently.
- 10.6 Using the table below as a guide, repeat for tubes #4 through #6.
- 10.7 Standard B<sub>0</sub> contains no protein and is the blank control.

Standard #	Sample to Dilute	Volume to Dilute ( $\mu$ L)	Volume of Diluent ( $\mu$ L)	Starting Conc. (ng/mL)	Final Conc. (pg/mL)
1	<b>Step 10.3</b>				<b>1,000,000</b>
2	<b>Standard #1</b>	100	900	1,000,000	100,000
3	<b>Standard #2</b>	100	900	100,000	10,000
4	<b>Standard #3</b>	100	900	10,000	1,000
5	<b>Standard #4</b>	100	900	1,000	100
6	<b>Standard #5</b>	100	900	100	10
B <sub>0</sub>	<b>None</b>		900	-	-



### 11. SAMPLE COLLECTION AND STORAGE

There are different protocols for isolating and purifying 14,15-DHET depending on the medium in which it is in. For optimal results follow the appropriate protocol based on the biological sample present.

#### 11.1 **14,15-DHET measurement in cells**

- 11.1.1 Collect and homogenize and/or sonicate the cells using a solution containing a final concentration of ~0.1 mM TPP (triphenylphosphine, 0.03-0.05 mg/mL). TPP is an antioxidant, which looks like a precipitate in samples because it does not easily dissolve. Before using the stored samples containing TPP, spin samples to separate the precipitated TPP from sample solution.
- 11.1.2 Acidify the whole homogenized cells with acetic acid to a pH of approximately 3-4. Measure using standard pH paper.
- 11.1.3 Extraction with ethyl acetate. Add an equal volume of ethyl acetate to the homogenized cells and vortex very well. Place the upper organic phase into a fresh clean tube after centrifugation. Then add another equal volume of ethyl acetate to the homogenized cells to start the second-time

extraction. It is strongly recommended that extraction is performed three times.

- 11.1.4 Evaporate the pooled ethyl acetate from the extractions until all has dried up under argon or nitrogen gas.
  - 11.1.5 **Saponification if needed (see below)**
  - 11.1.6 Add 10  $\mu\text{L}$  to 20  $\mu\text{L}$  ethanol, or N, N-dimethylformamide (DMF), to reconstitute the dried-down residue from above step. Add 500  $\mu\text{L}$  of 1X Sample Dilution Buffer (provided in kit). Load 100  $\mu\text{L}$  in each well, in triplicates, on the ELISA plate. (Note: We recommend measuring a different dilution of sample in attempt to fit the results to the standard curve. e.g., add 3 wells with 50  $\mu\text{L}$  of the rest of sample plus 50  $\mu\text{L}$  1x Sample Dilution Buffer, and 3 wells with 10  $\mu\text{L}$  of the rest of sample plus 90  $\mu\text{L}$  of 1x Sample Dilution Buffer.)
  - 11.1.7 Perform the ELISA for 14,15-DHET
- 11.2 **Saponification** (to cleave fatty acid from glycerol backbone):
- 11.2.1 Dissolve dried fatty acids (obtained from 3X ethyl acetate extractions) in 2 mL of 20% KOH solution (make working solution: 1 mL of 2 M KOH + 4 mL methanol so that the final conc. of KOH = 0.4 N).
  - 11.2.2 Vortex and incubate for 1 h at 50°C.
  - 11.2.3 Add 1.5 X H<sub>2</sub>O to the solution and adjust pH with 20% formic acid to pH~5.
  - 11.2.4 Re-extract the solution with ethyl acetate (1 part aqueous solution + 1 part ethyl acetate) and dry.
- 11.3 **14,15-DHET measurement in tissues**

- 11.3.1 Homogenize 1 g of tissue, 4 mL of H<sub>2</sub>O, and 0.01 mg TPP.
  - 11.3.2 Acidify the homogenate by adding 8 µL of acetic acid to each homogenate.
  - 11.3.3 Extract with an equal amount of ethyl acetate, vortex thoroughly, spin down, and collect the organic phase. Repeat this extraction twice more and combine all of the organic phases.
  - 11.3.4 Dry the organic phase with argon or nitrogen gas.
  - 11.3.5 **Saponification if needed (see Section 11.2)**
  - 11.3.6 Dissolve the dried residue from above step with ethanol or DMF. (Add approximately 20 µL of ethanol or DMF to reconstitute the dried-up residue.)
  - 11.3.7 Dilute further with 1X Sample Dilution Buffer: Add approximately 500 µL of 1x Sample Dilution Buffer and centrifuge at 10,000 rpm for five minutes at room temperature. The supernatant will be used for ELISA.
  - 11.3.8 Perform the ELISA for 14,15-DHET (according to the instructions of the manufacturer).
- 11.4 **14,15-DHET protocol for 1.0 mL plasma or serum**
- 11.4.1 Combine 1 mL of plasma (adjusted with approximately 12 µL of acetic acid to pH 4) and 1 mL of ethyl acetate. Vortex thoroughly. Centrifuge at 2,000 rpm for ten minutes at 22°C. Three phases should result:
    - 11.4.1.1 Upper organic phase – ethyl acetate phase (lipoproteins)
    - 11.4.1.2 Interphase – proteins
    - 11.4.1.3 Lower phase – aqueous phase
  - 11.4.2 Collect the upper organic phase (a) and set aside.

## ASSAY PREPARATION

- 11.4.3 Discard the interphase. Transfer the lower phase with a glass pipette to a new tube, and repeat the ethyl acetate extraction step 2 more times.
- 11.4.4 Evaporation of pooled organic phase: There should be approximately 3 mL of the ethyl acetate phase (a). Dry the pooled organic phase in a Speedvac to get the extracted sediment (b).
- 11.4.5 Saponification (to cleave fatty acid from glycerol backbone): Dissolve the dried residues (b) in 2 mL of 20% KOH solution (for preparation see 14,15-DHET measurement in cells). Vortex thoroughly and incubate for 1 h at 50°C. This will yield an aqueous solution (c).
- 11.4.6 Dilute 2 mL of the aqueous solution (c) with 3 mL of H<sub>2</sub>O. Adjust the pH using 20% formic acid (132 µL) to pH~5.5. Add ethyl acetate (1 part aqueous solution (c) + 1 part ethyl acetate), vortex thoroughly, and centrifuge at 2,000 rpm for ten minutes at 22°C. Repeat the procedure twice more using an equal volume of ethyl acetate per sample. Collect the upper phase containing saponified lipids.
- 11.4.7 Dry the pooled ethyl acetate upper phase (d) and dry in a Speedvac, yielding the dried sample-sediment (e). Store the sediment (e) at -20°C. For ELISA assay, dissolve the sediment (e) in 20 µL of ethanol, then add 380 µL of 1X Sample Dilution Buffer, pH 7.4. *(Please note that the 10X Sample Dilution Buffer that is supplied with the ELISA kit must be diluted 10-fold).*
- 11.4.8 When calculating the concentration, consider the dilution factor. In this case, 400 µL total sample volume from 1 mL plasma (2.5-fold concentration) you must divide your calculated result by 2.5.
- 11.4.9 Perform the ELISA for 14,15-DHET

### 11.5 **14,15-DHET protocol for 1.8 mL plasma or serum**

***(Following this procedure, the user will have approximately 70  $\mu$ L of material left-over from step 11.5.8. This material can be stored at -20°C and be used for a second measurement following a 5X dilution).***

- 10.5.1 Combine 1.8 mL of plasma (adjusted with approximately 20  $\mu$ L of acetic acid to pH 4) and 1.8 mL of ethyl acetate. Vortex thoroughly. Centrifuge at 2,000 rpm for ten minutes at 22°C. Three phases should result:
  - 11.5.1.1 Upper organic phase – ethyl acetate phase (lipoproteins)
  - 11.5.1.2 Interphase – proteins
  - 11.5.1.3 Lower phase – aqueous phase
- 11.5.2 Collect the upper organic phase (a) and set aside.
- 11.5.3 Discard the interphase. Transfer the lower phase with a glass pipette to a new tube, and repeat the ethyl acetate extraction step 2 more times.
- 11.5.4 Evaporation of pooled organic phase: There should be approximately 5-6 mL of the ethyl acetate phase (a). Dry the pooled organic phase in a Speedvac to get the extracted sediment (b).
- 11.5.5 Saponification (to cleave fatty acid from glycerol backbone): Dissolve the dried residues (b) in 2 mL of 20% KOH solution (for preparation see 14,15-DHET measurement in cells). Vortex thoroughly and incubate for 1 h at 50°C. This will yield an aqueous solution (c).
- 11.5.6 Dilute 2 mL of the aqueous solution (c) with 3 mL of H<sub>2</sub>O. Adjust the pH using 20% formic acid (132  $\mu$ L) to pH~5.5. Add ethyl acetate (1 part aqueous solution (c) + 1 part ethyl acetate), vortex thoroughly, and centrifuge at 2,000 rpm for ten minutes at 22°C. Repeat the procedure twice more

- using an equal volume of ethyl acetate per sample. Collect the upper phase containing saponified lipids.
- 11.5.7 Dry the pooled ethyl acetate upper phase (d) and dry in a Speedvac, yielding the dried sample-sediment (e). Store the sediment (e) at  $-20^{\circ}\text{C}$ . For ELISA assay, dissolve the sediment (e) in  $20\ \mu\text{L}$  of ethanol, then add  $130\ \mu\text{L}$  of 1X Sample Dilution Buffer.
  - 11.5.8 For the competitive 14,15-DHET ELISA, the above  $150\ \mu\text{L}$  sample needs to be further diluted: Dilute 1:4 (e.g.,  $80\ \mu\text{L}$  sample +  $320\ \mu\text{L}$  1x Sample Dilution Buffer). Check the final pH (should be pH 7.4). When calculating the concentration, consider the dilution factor. In this case,  $150\ \mu\text{L}$  total sample volume from  $1.8\ \text{mL}$  plasma (12-fold concentration) and then,  $80\ \mu\text{L}$  sample in  $400\ \mu\text{L}$  SDB (5-fold dilution). Since, the samples are concentrated 2.4-fold; to get the actual concentration, you must divide by 2.4.
  - 11.5.9 Perform the ELISA for 14,15-DHET (according to the instructions of the manufacturer).
- 11.6 14,15-DHET measurement in urine**
- 11.6.1 Extraction using ethyl acetate is not necessary. It is recommended that the urine sample be diluted 4-fold with 1X Sample Dilution Buffer and  $100\ \mu\text{L}$  of sample added directly to the ELISA plate well.
  - 11.6.2 However, extraction of urine with ethyl acetate can be performed if desired. See protocols for extraction with ethyl acetate above.
  - 11.6.3 When calculating the concentration, consider the dilution factor.

## 12. PLATE PREPARATION

- The 96 well plate included with this kit are supplied ready to use. It is not necessary to rinse the plate prior to adding reagents
- Unused well strips should be returned to the plate packet and stored at 4°C
- For each assay performed, a minimum of 2 wells must be used as a blank, omitting sample and conjugate from well addition. Another 2 wells must be used for a maximum binding control
- For statistical reasons, we recommend each standard and sample should be assayed with a minimum of two replicates (duplicates)

## **13. ASSAY PROCEDURE**

- **Equilibrate all materials and prepared reagents to room temperature prior to use**
  - **Please read the test protocol carefully before performing the assay. Result reliability depends on strict adherence to the test protocol as described**
  - **If performing the test on an automatic ELISA system we recommend increasing the washing steps from three to five and the volume of 1X Wash Buffer from 300  $\mu$ L to 350  $\mu$ L to avoid washing effects**
  - **Assay all standards, controls and samples in duplicate**
- 
- 13.1 Add 200  $\mu$ L of 1X Sample Dilution Buffer into the blank wells and 100  $\mu$ L of 1X Sample Dilution Buffer into maximum binding control wells.
  - 13.2 Add 100  $\mu$ L of each of the standards or samples into the appropriate wells.
  - 13.3 Add 100  $\mu$ L of the 1X-HRP conjugate in the all wells except the blank control wells.
  - 13.4 Incubate the plate at room temperature for two hours.
  - 13.5 Wash the plate three times with 400  $\mu$ L of 1X Wash Buffer per well.
  - 13.6 After the last of the three wash cycles pat the inverted plate dry onto some paper towels.
  - 13.7 Add 200  $\mu$ L of the TMB substrate to all of the wells.
  - 13.8 Incubate the plate at room temperature for 15-30 minutes.
  - 13.9 Add 50  $\mu$ L of 2 N sulfuric acid to all of the wells.
  - 13.10 Read the plate at 450 nm.


## 14. CALCULATIONS

If data redaction software is not available on your plate reader then the results can be obtained manually as follows:

- 14.1 Average the absorbance (Abs) readings from the blank wells and subtract that value from each well of the plate to obtain the corrected readings. (Note: Some plate readers do this automatically. Consult the user manual of your plate reader.)
- 14.2 Average the corrected absorbance readings from the maximum binding control wells. This is your maximum binding.
- 14.3 Calculate the % Abs for Standard 1 by averaging the corrected absorbance of the two wells; divide the average by the Maximum Binding Control well average absorbance, then multiply by 100. Repeat this formula for the remaining standards.
- 14.4 Plot the % Abs versus the concentration of 14,15-DHET from the standards using semi-log paper.
- 14.5 Calculate the % Abs for the samples and determine the concentrations, utilizing the standard curve.
- 14.6 Multiply the concentrations obtained for each of the samples by their corresponding dilution factor.

## 15. TYPICAL DATA

**TYPICAL STANDARD CURVE** – Data provided for demonstration purposes only. A new standard curve must be generated for each assay performed. The data shown here is an example of typical results obtained using the Abcam's 14,15-DHET ELISA kit. These results are only a guideline, and should not be used to determine values from your samples.



Conc. (pg/mL)	% (B/B <sub>0</sub> )
10	85.4
100	68.2
1,000	47.3
10,000	27.9
100,000	11.0
1,000,000	5.4

## 16. ASSAY SPECIFICITY

The specificity of the 14,15-DHET ELISA was investigated using authentic 14,15-DHET and a panel of eicosanoids.

Eicosanoid	Reactivity
14,15-DHET	100.00%
8,9-DHET	3.30%
11,12-DHET	3.30%
14,15-EET	1.5 %
15(s) HETE	1.00%
8,9-EET	0.40%
5(s)15(s)DiHETE	0.20%
11,12-EET	0.05%
Arachidonic Acid	0.05%
5,6-DHET	0.02%
5,6-EET	0.02%
Thromboxane B <sub>2</sub>	0.02%
PGE <sub>2</sub>	<0.01 %
PGF <sub>2a</sub>	<0.01 %
6-keto-PGF <sub>1a</sub>	<0.01 %

## SENSITIVITY-

The calculated minimal detectable (MDD) dose is 3 pg/mL. The MDD was determined by calculating the mean of zero standard replicates.

## 17. TROUBLESHOOTING

<b>Problem</b>	<b>Cause</b>	<b>Solution</b>
No color present in standard wells	The HRP conjugate was not added	Redo the assay and add the conjugate at the proper step
	The HRP conjugate was not incubated for the proper time	Redo the assay and incubate for the proper time
No color in any wells	The TMB substrate was not added	Add substrate
	The TMB conjugate was not incubated for the proper time	Continue incubation until desired color is reached

## RESOURCES

<b>Problem</b>	<b>Cause</b>	<b>Solution</b>
The color is faint	One or all of the incubation times were cut short	Redo the assay with the proper incubation times
	The TMB substrate was not warmed up to room temperature	Redo the assay making sure all reagents are at room temperature
	The lab is too cold	Be sure the lab temperature is between 21-27°C and redo the assay
The background color is very high	The TMB substrate has been contaminated	Redo the assay with a fresh bottle of substrate
Scattered OD obtained from sample	Incorrect loading of samples	Redo assay using an 8 channel pipetman making sure the 8 channels are equal volume while loading

18. NOTES

**UK, EU and ROW**

Email: [technical@abcam.com](mailto:technical@abcam.com) | Tel: +44-(0)1223-696000

**Austria**

Email: [wissenschaftlicherdienst@abcam.com](mailto:wissenschaftlicherdienst@abcam.com) | Tel: 019-288-259

**France**

Email: [supportscientifique@abcam.com](mailto:supportscientifique@abcam.com) | Tel: 01-46-94-62-96

**Germany**

Email: [wissenschaftlicherdienst@abcam.com](mailto:wissenschaftlicherdienst@abcam.com) | Tel: 030-896-779-154

**Spain**

Email: [soportecientifico@abcam.com](mailto:soportecientifico@abcam.com) | Tel: 911-146-554

**Switzerland**

Email: [technical@abcam.com](mailto:technical@abcam.com)

Tel (Deutsch): 0435-016-424 | Tel (Français): 0615-000-530

**US and Latin America**

Email: [us.technical@abcam.com](mailto:us.technical@abcam.com) | Tel: 888-77-ABCAM (22226)

**Canada**

Email: [ca.technical@abcam.com](mailto:ca.technical@abcam.com) | Tel: 877-749-8807

**China and Asia Pacific**

Email: [hk.technical@abcam.com](mailto:hk.technical@abcam.com) | Tel: 108008523689 (中國聯通)

**Japan**

Email: [technical@abcam.co.jp](mailto:technical@abcam.co.jp) | Tel: +81-(0)3-6231-0940

[www.abcam.com](http://www.abcam.com) | [www.abcam.cn](http://www.abcam.cn) | [www.abcam.co.jp](http://www.abcam.co.jp)