

Version 1 Last updated 6 December 2018

ab242287

Alkaline Phosphatase Staining Kit (Purple)

[View Alkaline Phosphatase Staining Kit datasheet:
www.abcam.com/ab242287](https://www.abcam.com/ab242287)

For the for monitoring of ES cell undifferentiation/ differentiation.

This product is for research use only and is not intended for diagnostic use.

Table of Contents

1. Overview	3
2. Protocol Summary	4
3. General guidelines, precautions, and troubleshooting	5
4. Materials Supplied, and Storage and Stability	6
5. Materials Required, Not Supplied	6
6. Reagent Preparation	7
7. Assay Procedure	8
8. Typical Data	9
9. Notes	10

1. Overview

Like mouse Embryonic stem (ES) cells, human ES cells express high levels of membrane alkaline phosphatase (AP) and Oct-4, a transcriptional factor critical to inner cell mass (ICM) and germline formation. However, unlike mouse ES cells, hES cells do not express stage-specific embryonic antigen (SSEA-1). In addition, prolonged propagation of hES cells is typically achieved by co-culture with primary mouse embryonic fibroblasts (MEFs) serving as feeder cells. Human ES cell lines are not able to maintain their undifferentiated state in the absence of supporting feeder layer cells, even when exogenous cytokines such as leukemia inhibitory factor (LIF) and gelatin-coated plates are used.

Marker Name	Mouse ES Cells	Mouse EG Cells	Human ES Cells	Human EG Cells	Human EC Cells
AP	Yes	Yes	Yes	Yes	Yes
SSEA-1	Yes	Yes	No	Yes	No
SSEA-4	No	No	Yes	Yes	Yes
TRA-1-60	No	No	Yes	Yes	Yes
TRA-1-81	No	No	Yes	Yes	Yes
Oct-4	yes	yes	yes	unknown	Yes

ES Cell = Embryonic stem cell

EG Cell = Embryonic germ cell

EC Cell = Embryonic carcinoma cell

Figure 1: Comparison of mouse and human pluripotent stem cells.

Although stem cells from different origins require different growth conditions for self-renewal and display different cell surface markers, AP is the most widely used stem cell marker.

Alkaline Phosphatase Staining Kit (Purple) (ab242287) provides an efficient system for monitoring ES cell undifferentiation/ differentiation through AP activity by immunocytochemistry staining.

2. Protocol Summary

Culture mouse ES cells in medium containing LIF; human ES cells on a MEF feeder layer.



Aspirate medium and wash with PBST.



Aspirate PBST and add Fixing Solution, incubate at room temperature for 2 mins.



Remove Fixing Solution and wash with PBST.



Aspirate wash and add AP Staining Solution, incubate at room temperature for 15-30 mins, protected from light.



Remove AP Staining Solution and wash with 1X PBS.



Either overlay cells with 1X PBS/ 20% glycerol at store at +4°C or count the purple stained colonies vs. colorless colonies using a light microscope.

3. General guidelines, precautions, and troubleshooting

- Please observe safe laboratory practice and consult the safety datasheet.
- For general guidelines, precautions, limitations on the use of our assay kits and general assay troubleshooting tips, particularly for first time users, please consult our guide:
www.abcam.com/assaykitguidelines
- For typical data produced using the assay, please see the assay kit datasheet on our website.

4. Materials Supplied, and Storage and Stability

- Store kit at 4°C immediately upon receipt and check below for storage for individual components. Kit can be stored for 1 year from receipt, if components have not been reconstituted.
- Aliquot components in working volumes before storing at the recommended temperature.
- Avoid repeated freeze-thaws of reagents.

Item	Quantity	Storage condition
Fixing Solution	50 mL	4°C
AP Staining Solution	40 mL	4°C

5. Materials Required, Not Supplied

These materials are not included in the kit, but will be required to successfully perform this assay:

- Human or mouse embryonic stem cells and culture medium
- 1X PBS
- 1X PBST (1X PBS containing 0.05% Tween-20)
- Deionized water
- Light microscope

6. Reagent Preparation

- Equilibrate all reagents to room temperature (18-25°C) prior to use. Before using the kit, spin tubes and bring down all components to the bottom of tubes.
- Prepare only as much reagent as is needed on the day of the experiment.
- Any components not listed here are ready to use as supplied.

7. Assay Procedure

7.1 24-Well Plate:

- 7.1.1 Culture mouse ES cells in medium containing LIF; alternatively, culture human ES cells on a MEF feeder layer.
- 7.1.2 Gently aspirate the medium from the ES cells and wash the cells with 1 mL of 1X PBST. Aspirate the wash solution.
- 7.1.3 Add Fixing Solution to the cells, 0.4 mL per well for a 24-well plate. Incubate at room temperature for 2 mins.
- 7.1.4 Remove the fixing solution and wash the fixed cells twice with 1 mL of 1X PBST.
- 7.1.5 Aspirate the final wash and add 0.4 mL per well AP Staining Solution.
- 7.1.6 Incubate the cells at room temperature for 15-30 mins, protected from light.
- 7.1.7 Remove the AP Staining Solution, and then wash the stained cells twice with 1 mL of 1X PBS. Store cells in 1X PBS at 4°C. For long-term storage, overlay the cells with 1X PBS containing 20% Glycerol. Store at 4°C.
- 7.1.8 Count the purple stained cell colonies (undifferentiated ES cells) vs. colorless colonies (differentiated ES cells) using a light microscope.

8. Typical Data

The following figure demonstrates typical results with AP Staining Solution. One should use the data below for reference only. This data should not be used to interpret actual results.




Figure 1. AP staining of ES Cells. Murine embryonic stem cells (ES-D3) are maintained in an undifferentiated stage on gelatin-coated dishes in the presence of LIF. AP staining of undifferentiated ES cells was performed as described in the Assay Protocol.

9. Notes

Technical Support

Copyright © 2018 Abcam. All Rights Reserved. The Abcam logo is a registered trademark. All information / detail is correct at time of going to print.

Austria

wissenschaftlicherdienst@abcam.com | 019-288-259

France

supportscientifique@abcam.com | 01.46.94.62.96

Germany

wissenschaftlicherdienst@abcam.com | 030-896-779-154

Spain

soportecientifico@abcam.com | 91-114-65-60

Switzerland

technical@abcam.com

Deutsch: 043-501-64-24 | Français: 061-500-05-30

UK, EU and ROW

technical@abcam.com | +44(0)1223-696000

Canada

ca.technical@abcam.com | 877-749-8807

US and Latin America

us.technical@abcam.com | 888-772-2226

Asia Pacific

hk.technical@abcam.com | (852) 2603-6823

China

cn.technical@abcam.com | +86 21 2070 0500 | 400 921 0189

Japan

technical@abcam.co.jp | +81-(0)3-6231-0940

Singapore

sg.technical@abcam.com | 800 188-5244

Australia

au.technical@abcam.com | +61-(0)3-8652-1450

New Zealand

nz.technical@abc.com | +64-(0)9-909-7829