

ab325370 – Human FGL1 SimpleStep ELISA® Kit, Chemiluminescent

For the quantitative measurement of FGL1 in human serum, plasma (heparin), plasma (EDTA), plasma (citrate), and cell culture supernatant.

For research use only - not intended for diagnostic use.

Storage and Stability: Store kit at 2-8°C immediately upon receipt. Refer to list of materials supplied for storage conditions of individual components. Observe the storage conditions for individual prepared components in the Standard Preparation and Reagent preparation sections.

Limitations: Typical Standard Curve, Sensitivity, and individual donor serum were determined using Chemiluminescent detection. Unless stated, all remaining data was collected using the colorimetric version of this kit (ab284622).

Materials Supplied

Item	Quantity 1 x 96 well	Quantity 10 x 96 well	Storage Condition
Human FGL1 Capture Antibody 10X	600 µL	10 x 600 µL	+4°C
Human FGL1 Detector Antibody 10X	600 µL	10 x 600 µL	+4°C
Human FGL1 Lyophilized Recombinant Protein	2 Vials	10 x 2 Vials	+4°C
Antibody Diluent 5BI	6 mL	10 x 6 mL	+4°C
Sample Diluent NS	50 mL	2 x 250 mL	+4°C
Wash Buffer PT 10X	20 mL	200 mL	+4°C
ChemiHRP Reagent A	3 mL	50 mL	+4°C
ChemiHRP Reagent B	3 mL	50 mL	+4°C
SimpleStep Pre-Coated Black 96-Well Microplate	96 Wells	10 x 96 Wells	+4°C
Plate Seal	1	10	+4°C

Materials Required, Not Supplied

These materials are not included in the kit, but will be required to successfully utilize this assay:

Luminometer with the following settings: 0.5-1 second/well read time; summation mode (all wavelengths)

Deionized water.

Multi- and single-channel pipettes.

Tubes for standard dilution.

Plate shaker for all incubation steps.

Optional: Phenylmethylsulfonyl Fluoride (PMSF) (or other protease inhibitors).

Reagent Preparation

Equilibrate all reagents to room temperature (18-25°C) prior to use. The kit contains enough reagents for 96 wells. The sample volumes below are sufficient for 48 wells (6 x 8-well strips); adjust volumes as needed for the number of strips in your experiment.

Prepare only as much reagent as is needed on the day of the experiment. Capture and Detector Antibodies have only been tested for stability in the provided 10X formulations.

1X Wash Buffer PT: Prepare 1X Wash Buffer PT by diluting Wash Buffer PT 10X with deionized water. To make 50 mL 1X Wash Buffer PT combine 5 mL Wash Buffer PT 10X with 45 mL deionized water. Mix thoroughly and gently.

Antibody Cocktail: Prepare Antibody Cocktail by diluting the capture and detector antibodies in Antibody Diluent 5BI. To make 3 mL of the Antibody Cocktail combine 300 µL 10X Capture Antibody and 300 µL 10X Detector Antibody with 2.4 mL Antibody Diluent 5BI. Mix thoroughly and gently.

Lumi HRP Development Solution: Just prior to use, prepare Lumi HRP Development Solution by mixing equal volume of the ChemiHRP Reagent A and the ChemiHRP Reagent B. To make 3 mL of the Lumi HRP Development Solution combine 1.5 mL of ChemiHRP Reagent A and 1.5 mL of ChemiHRP Reagent B. Mix thoroughly and gently by inversion or slow pipetting (Avoid shaking or vortexing). Protect the prepared solution from light until use.

Standard Preparation

Always prepare a fresh set of standards for every use. Discard working standard dilutions after use as they do not store well. The following section describes the preparation of a standard curve for duplicate measurements (recommended).

1. Reconstitute the FGL1 standard sample by adding the volume of Sample Diluent NS indicated on the protein vial label. Hold at room temperature for 10 minutes. Mix thoroughly and gently. This is the 10,000 pg/mL **Stock Standard** Solution.
2. Label nine tubes, Standards 1– 9.
3. Add 150 µL of Sample Diluent NS into numbers 2-9.
4. Use the **Stock Standard** to prepare the following dilution series. Standard #9 contains no protein and is the Blank control:

Standard #	Dilution Sample	Volume to Dilute (µL)	Volume of Diluent (µL)	Starting Conc. (pg/mL)	Final Conc. (pg/mL)
1	Stock Standard	225	0	10,000	10,000
2	Standard#1	75	150	10,000	3,333.3
3	Standard#2	75	150	3,333.3	1,111.1
4	Standard#3	75	150	1,111.1	370.4
5	Standard#4	75	150	370.4	123.5
6	Standard#5	75	150	123.5	41.2
7	Standard#6	75	150	41.2	13.7
8	Standard#7	75	150	13.7	4.6
9	Blank Control	0	150	0	0

Sample Preparation

Typical Sample Dynamic Range	
Sample Type	Range
Serum	0.063 – 1%
Plasma – Citrate	0.016 – 0.5%
Plasma – EDTA	0.016 – 0.5%
Plasma – Heparin	0.016 – 0.5%
HepG2 Cell Culture Supernatant	0.031 – 0.5%

Serum Samples should be collected into a serum separator tube. After clot formation, centrifuge samples at 2,000 x g for 10 minutes and collect serum. Dilute samples at least 1:100 into Sample Diluent NS and assay. Store un-diluted serum at -20°C or below. Avoid repeated freeze-thaw cycles.

Plasma Collect plasma using citrate, EDTA, or heparin. Centrifuge samples at 2,000 x g for 10 minutes. Dilute samples at least 1:200 into Sample Diluent NS and assay. Store un-diluted plasma samples at -20°C or below for up to 3 months. Avoid repeated freeze-thaw cycles.

Cell Culture Supernatants Centrifuge cell culture media at 2,000 x g for 10 minutes to remove debris. Collect supernatants and assay. Or dilute samples at least 1:200 into Sample Diluent NS and assay. Store un-diluted samples at -20°C or below. Avoid repeated freeze-thaw cycles.

Plate Preparation

The 96 well plate strips included with this kit are supplied ready to use. It is not necessary to rinse the plate prior to adding reagents.

Unused plate strips should be immediately returned to the foil pouch containing the desiccant pack, resealed and stored at 4°C.

For each assay performed, a minimum of two wells must be used as the zero control.

For statistical reasons, we recommend each sample should be assayed with a minimum of two replicates (duplicates).

Differences in well absorbance or “edge effects” have not been observed with this assay.

Assay Procedure

Equilibrate all materials and prepared reagents to room temperature prior to use.

We recommend that you assay all standards, controls and samples in duplicate.

1. Prepare all reagents, working standards, and samples as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, reseal and return to 4°C storage.

3. Add 50 µL of all sample or standard to appropriate wells.
4. Add 50 µL of the Antibody Cocktail to each well.
5. Seal the plate and incubate for 30 minutes at room temperature on a plate shaker set to 750 rpm.
6. Wash each well with 3 x 350 µL 1X Wash Buffer PT. Wash by aspirating or decanting from wells then dispensing 350 µL 1X Wash Buffer PT into each well. Wash Buffer PT should remain in wells for at least 30 seconds. Complete removal of liquid at each step is essential for good performance. After the last wash invert the plate and tap gently against clean paper towels to remove excess liquid.
7. Add 50 µL of prepared Lumi HRP Development Solution to each well and incubate for 1 minute in the dark on a plate shaker set to 750 rpm. Further optimization of incubation time vs signal strength can be performed if needed. Avoid introducing bubbles into the wells.
8. Measure the produced light of each well using a microplate luminometer with the following settings: 0.5-1 second/well read time in summation mode (all wavelengths). Relative light unit (RLU) readings may vary between luminometer models. It is recommended to configure instrument settings according to the manufacturer's specifications. Note: Relative light unit (RLU) values may change over the course of the 15-minute reading window.
9. Analyze the data as described below.

Mode:	Luminescence
Instrument settings:	Endpoint
Detection Mode:	All wavelengths
Read Time:	0.5-1 sec
Read:	Top

Note For microplate readers with Pre-Read Optimization option, the Read Height as well as Microplate Optimization is recommended before the first read.

Download our ELISA guide for technical hints, results, calculation, and troubleshooting tips:

www.abcam.com/protocols/the-complete-elisa-guide

Technical Support

Copyright © 2026 Abcam. All Rights Reserved. The Abcam logo is a registered trademark. All information / detail is correct at time of going to print.

For all technical or commercial enquiries please go to:

<https://www.abcam.com/en-us/contact-us>

<https://www.abcam.cn/contact-us> (China)

<https://www.abcam.co.jp/contact-us> (Japan)

ab325370 – Human FGL1 SimpleStep ELISA® Kit, Chemiluminescent

Additional information

ASSAY SPECIFICITY

This kit is designed for the quantification of human FGL1.

The standard protein in this kit is mature full length human FGL1.

Native signal was detected in serum, plasma (heparin), plasma (EDTA), plasma (citrate), and cell culture supernatant.

Saliva, urine, milk, and CSF samples have not been tested with this kit.

This kit is incompatible with cell extract and tissue extract samples.

For the measurement of FGL1 in cell extract and tissue extract sample types, use colorimetric Human FGL1 ELISA kit ab284622.

SPECIES REACTIVITY

Other species reactivity was determined by measuring 0.25% serum samples of various species, interpolating the protein concentrations from the human standard curve, and expressing the interpolated concentrations as a percentage of the protein concentration in human serum assayed at the same dilution.

Reactivity was determined for the following species: Mouse (57%), Rat (15%).

Reactivity < 3% was determined for the following species: Monkey, Cow

Other species reactivity not determined.

CALCULATION

- Preconfigured protocols are available when using SoftMax Pro software from Molecular Devices.
- Calculate the average chemiluminescence value for the blank control (zero) standards. Subtract the average blank control standard chemiluminescence value from all other chemiluminescence values.
- Create a standard curve by plotting the average blank control subtracted chemiluminescence value for each standard concentration (y-axis) against the target protein concentration (x-axis) of the standard. Use graphing software to draw the best smooth curve through these points to construct the standard curve.
- Note: Most chemiluminescence reader software or graphing software will plot these values and fit a curve to the data. A four-parameter curve fit (4PL) is often the best choice; however, other algorithms (e.g. linear, semi-log, log/log, 4-parameter logistic) can also be tested to determine if it provides a better curve fit to the standard values.
- Determine the concentration of the target protein in the sample by interpolating the blank control subtracted chemiluminescence values against the standard curve. Multiply the resulting value by the appropriate sample dilution factor, if used, to obtain the concentration of target protein in the sample.

- Samples generating chemiluminescence values greater than that of the highest standard should be further diluted and reanalyzed. Similarly, samples which measure at chemiluminescence values less than that of the lowest standard should be retested in a less dilute form.

TYPICAL DATA

Typical standard curve – data provided for demonstration purposes only. A new standard curve must be generated for each assay performed

Standard Curve Measurements			
Concentration (pg/mL)	RFU		Mean RFU
	1	2	
0	809	788	799
4.6	3,008	2,741	2,875
13.7	6,561	6,157	6,359
41.2	15,080	17,777	16,429
123.5	52,420	51,287	51,854
370.4	150,370	161,560	155,965
1,111.1	480,930	467,050	473,990
3,333.3	1,431,800	1,431,700	1,431,750
10,000	4,155,700	4,104,000	4,129,850

Table 1. Example of human FGL1 standard curve in Sample Diluent NS. The FGL1 standard curve was prepared as described in the Standard Preparation section. The table shows raw data values.

TYPICAL SAMPLE VALUES

Sensitivity:

The calculated minimal detectable dose (MDD) is 1.7 pg/mL. The MDD was determined by calculating the mean of zero standard replicates (n=8) and adding 2 standard deviations then extrapolating the corresponding concentration.

Recovery

Three concentrations of FGL1 were spiked in duplicate to the indicated biological matrix to evaluate signal recovery in the working range of the assay.

Sample Type	Average % Recovery	Range (%)
0.25% Serum	106	95 - 117
0.25% Plasma – Citrate	96	87 - 103
0.25% Plasma – EDTA	108	102 - 113
0.25% Plasma – Heparin	95	92 - 99
0.25% HepG2 Supernatant	105	96 - 118

Linearity of Dilution

Linearity of dilution is determined based on interpolated values from the standard curve. Linearity of dilution defines a sample concentration interval in which interpolated target concentrations are directly proportional to sample dilution.

Native FGL1 was measured in the following biological samples in a 2-fold dilution series. Sample dilutions are made in Sample Diluent NS.

Dilution Factor	Interpolated value	1% Human Serum	0.5% Human Plasma (Citrate)	0.5% Human Plasma (EDTA)	0.5% Human Plasma (Heparin)	0.5% HepG2 Cell Supernatant
Undiluted	pg/mL	786.35	1,952.91	2,336.89	2,123.77	1,681.04
	% Expected value	100	100	100	100	100
2	pg/mL	409.97	1,008.67	1,147.30	1,069.21	775.00
	% Expected value	104	103	98	101	92
4	pg/mL	218.47	499.62	582.99	527.71	413.76
	% Expected value	111	102	100	99	98
8	pg/mL	111.37	261.70	305.45	264.77	206.48
	% Expected value	113	107	105	100	98
16	pg/mL	NL	130.99	145.74	139.22	105.07
	% Expected value	NL	107	100	105	100

Precision

Mean coefficient of variations of interpolated values of FGL1 from three concentrations of serum within the working range of the assay.

	Intra-assay	Inter-assay
N=	8	3
CV (%)	4.3	4.5

Download our ELISA guide for technical hints, results, calculation, and troubleshooting tips:

www.abcam.com/protocols/the-complete-elisa-guide

Technical Support

Copyright © 2026 Abcam. All Rights Reserved. The Abcam logo is a registered trademark. All information / detail is correct at time of going to print.

For all technical or commercial enquiries please go to:

<https://www.abcam.com/en-us/contact-us>

<https://www.abcam.cn/contact-us> (China)

<https://www.abcam.co.jp/contact-us> (Japan)

Copyright © 2026 Abcam, All Rights Reserved. All information / detail is correct at time of going to print.

Version 2 | 2026-08-07