


Anti-Histone H4 (acetyl K16) antibody [EPR1004] - ChIP Grade (ab109463)


Target	Histone H4 (acetyl K16)
Sample	HeLa
Treatment	None
Antibody amount	2ug
Cell number	2.5X10^5
Sequencing	NovaSeq 6000 (PE 150)
Peak calling	MACS 2.1.0
Reference genome	hg38

Chromosome 5




Copyright (c) 2021 Abcam plc.

Chromosome 5



Copyright (c) 2021 Abcam plc.

Chromosome 2



Copyright (c) 2021 Abcam plc.