



Anti-SP1 antibody [EPR6662(B)] - (ab124804)

Application	ChIC/CUT&RUN-seq
Target	SP1
Knock Out (KO) cell line	Human SP1 knockout HeLa cell line (ab265519)
Wild Type (WT) cell line	wild-type HeLa cell line (ab255928)
Treatment	None
Antibody amount	5µg
Cell number	2.5X10 ⁵
Antibody used for Assay Quality Control	Anti-CTCF(ab188408)
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	hg38
Reads	10 million

Chromosome 12 zoom out




ab124804 Chromosome 12




Copyright (c) 2022 Abcam plc.

ab124804 Chromosome 5




Copyright (c) 2022 Abcam plc.

ab124804 Chromosome 14



Copyright (c) 2022 Abcam plc.

ab124804 Chromosome 7



Copyright (c) 2022 Abcam plc.