



Anti-KDM1/LSD1 antibody [EPR6825] - ChIP Grade (ab129195)

Application	ChIC/CUT&RUN-seq
Target	KDM1/LSD1
Knock Out (KO) cell line	KDM1A (LSD1) knockout HeLa cell line (ab265790)
Wild Type (WT) cell line	wild-type HeLa cell line (ab255928)
Treatment	None
Antibody amount	5µg
Cell number	2.5X10 ⁵
Antibody used for Assay Quality Control	Anti-CTCF(ab188408)
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	hg38
Reads	10 million

Chromosome 9 zoom out




ab129195 Chromosome 9




Copyright (c) 2022 Abcam plc.

ab129195 Chromosome 1




Copyright (c) 2022 Abcam plc.

ab129195 Chromosome 1



Copyright (c) 2022 Abcam plc.

ab129195 Chromosome 10



Copyright (c) 2022 Abcam plc.