


Anti-Histone H3 (mono methyl K4) antibody [ERP16597] - ChIP Grade (ab176877)


|                  |                             |
|------------------|-----------------------------|
| Target           | Histone H3 (mono methyl K4) |
| Sample           | HeLa                        |
| Treatment        | None                        |
| Antibody amount  | 2ug                         |
| Cell number      | 2x 10^5                     |
| Sequencing       | NovaSeq 6000 (PE 150)       |
| Peak calling     | MACS 2.1.0                  |
| Reference genome | hg38                        |

Chromosome 7




Copyright (c) 2021 Abcam plc.

Chromosome 12




Copyright (c) 2021 Abcam plc.

Chromosome 12



Copyright (c) 2021 Abcam plc.

Chromosome 6



Copyright (c) 2021 Abcam plc.