

Peptide array: H3 K4me2 (ab176878)

Area under curve (AUC) data for all peptides on the pad. Antigen-containing peptides are highlighted in green.

Cell	AUC	Peptide
A1	8.5	ab10112 - H4 K91ac
A2	6.6	ab10121 - H2B K85ac
A3	9.5	ab10127 - H2A K95me2
A4	1.8	ab10138 - H1.4 T17p
A5	5.1	ab10139 - H1.4 T146p
A6	6.5	ab113997 - H2A.X Y142p
A7	7.7	ab11477 - H3 S10p
A8	6.6	ab12946 - H2B K5me3
A9	5.1	ab12949 - H2B K23me2
A10	6.2	ab12952 - H2A K127me1
A11	7.9	ab13211 - H2B K5me1
A12	10.7	ab1340 - H3 K4me1
A13	7.8	ab13570 - H2B K43me1
A14	7.8	ab13835 - HA
A15	5.8	ab13837 - Myc
A16	2.9	ab140408 - ScH2B T129p
A17	179.7	ab14103 - H3 T6p
A18	5.9	ab14724 - H4 S1p
A19	7.3	ab14799 - H3 T32p
A20	6.3	ab14943 - His
A21	7.6	ab14952 - H2A R77me2s
A22	5.4	ab14964 - H4 K20me2
A23	6.1	ab14997 - ScH2A K21ac
A24	6.5	ab154782 - H3 R26me2s
A25	5.4	ab15591 - H3 K14ac
B1	5.4	ab15644 - H3 K9me3 + S10p
B2	6.7	ab15645 - H2A.X S139p
B3	7.6	ab15824 - H4 K8ac
B4	4.3	ab15829 - V5
B5	9.1	ab16635 - H3 K9ac
B6	6.9	ab166685 - H2B S112GlcNAc
B7	7.1	ab166688 - H2B S36GlcNAc
B8	7.3	ab16935 - H3 R17me2a
B9	8.8	ab17043 - H4 K20me1
B10	7.6	ab17162 - H3 P16hyd + P30hyd
B11	8.4	ab17567 - H4 K20me3
B12	1.9	ab17587 - H1.4 K25me3
B13	9.7	ab175905 - H2B K11ac
B14	5.3	ab176211 - H2B K23ac
B15	2.8	ab17632 - H4 K5bio
B16	6.7	ab1771 - H3 K9me1
B17	8.4	ab1773 - H3 K9me3
B18	7.5	ab17770 - H4 R3me1
B19	5.7	ab1781 - H3 K27me2
B20	10.5	ab1783 - H3 K36me1
B21	17.1	ab1784 - H3 K36me2
B22	8.7	ab1785 - H3 K36me3
B23	11.0	ab178667 - ScH2A unmod(1-30)
B24	5.2	ab178668 - ScH2A unmod(102-132)
B25	11.0	ab178670 - H2A.Z unmod(1-30)
C1	6.9	ab178833 - Flag
C2	5.4	ab178835 - H2B R99mutO
C3	4.8	ab178836 - H4 R23mutO
C4	7.5	ab178837 - H4 G48mutD
C5	2.7	ab178987 - H1.4 unmod(1-30)
C6	1.8	ab178988 - H1.4 S1p
C7	1.9	ab178989 - H1.4 T3p
C8	1.7	ab178990 - H1.4 S1p + T3p
C9	1.4	ab178991 - H1.4 K25me1
C10	3.4	ab178992 - H1.4 K25me3 + S26p
C11	3.6	ab178993 - H1.4 unmod(130-160)

C12	12.3	ab178994 - H2A unmod(1-30)
C13	3.4	ab178995 - H2A S1p
C14	3.9	ab178996 - H2A S1p + K5ac
C15	6.8	ab178997 - H2A S1p + R3me2a + K5ac
C16	7.9	ab178998 - H2A S1p + R3me2s + K5ac
C17	8.0	ab178999 - H2A S1p + R3cit + K5ac
C18	7.8	ab179000 - H2A S1ac + S1p + K5ac + K9ac + K13ac + K15ac
C19	5.3	ab179001 - H2A S1p + K5ac + K9ac + K13ac + K15ac
C20	6.9	ab179002 - H2A R3cit
C21	7.1	ab179003 - H2A R3cit + K5ac
C22	7.9	ab179004 - H2A R3me2a + K5ac
C23	8.3	ab179005 - H2A R3me2s + K5ac
C24	10.6	ab179006 - H2A S1ac + K5ac + K9ac + K13ac + K15ac
C25	7.4	ab179007 - H2A K5ac + K9ac + K13ac + K15ac
D1	1.8	ab179008 - H2A K9me2
D2	0.8	ab179009 - H2A K9me3
D3	2.1	ab179010 - H2A R11me1
D4	1.5	ab179011 - H2A R11me2s
D5	2.7	ab179012 - H2A R11me2a
D6	1.7	ab179013 - H2A K13ac
D7	1.8	ab179014 - H2A K15ac
D8	2.2	ab179015 - H2A R17me1
D9	3.9	ab179016 - H2A R17me2s
D10	5.1	ab179017 - H2A R17me2a
D11	17.5	ab179018 - H2A unmod(20-50)
D12	4.8	ab179019 - H2A P26hyd
D13	6.2	ab179020 - H2A R29me1
D14	5.7	ab179021 - H2A R29me2s
D15	4.2	ab179022 - H2A R29me2a
D16	2.0	ab179023 - H2A K36ac
D17	6.3	ab179024 - H2A K36cr
D18	7.2	ab179025 - H2A Y39hyd
D19	7.9	ab179026 - H2A R42me1
D20	6.5	ab179027 - H2A R42me2s
D21	4.5	ab179028 - H2A R42me2a
D22	6.9	ab179029 - H2A unmod(71-103)
D23	8.2	ab179030 - H2A R77me1
D24	5.8	ab179031 - H2A R77me2a
D25	6.2	ab179032 - H2A R88me1
E1	6.8	ab179033 - H2A R88me2s
E2	8.8	ab179034 - H2A R88me2a
E3	5.6	ab179035 - H2A K95me1
E4	6.7	ab179036 - H2A K95me3
E5	10.4	ab179037 - H2A unmod(92-123)
E6	2.7	ab179038 - H2A T101GlcNAc
E7	2.0	ab179039 - H2A Q104me1
E8	2.2	ab179040 - H2A unmod(112-129)
E9	1.9	ab179041 - H2A K118me1
E10	2.4	ab179042 - H2A K118me2
E11	2.8	ab179043 - H2A K118me3
E12	2.4	ab179044 - H2A K118fo
E13	2.0	ab179045 - H2A K118cr
E14	4.0	ab179046 - H2A K119cr
E15	6.9	ab179047 - H2A K118ubiq
E16	2.4	ab179048 - H2A K119ubiq
E17	4.5	ab179049 - H2A T120p
E18	7.3	ab179050 - H2A K125me1
E19	4.3	ab179051 - H2A K125me2
E20	7.5	ab179052 - H2A K125cr
E21	21.1	ab179053 - H2A.X unmod(1-30)
E22	9.4	ab179054 - H2A.X K5ac
E23	8.1	ab179055 - H2A.X unmod(114-142)
E24	7.1	ab179056 - H2A.X K118ubiq
E25	10.0	ab179057 - H2A.X K119ubiq

F1	9.5	ab179058 - H2A.X T120p
F2	7.6	ab179059 - H2A.X K127ubiq
F3	8.2	ab179060 - H2A.X S139p + Y142p
F4	3.6	ab179061 - H2A.Z K4ac
F5	10.2	ab179062 - H2A.Z K7ac
F6	6.4	ab179063 - H2A.Z K11ac
F7	9.3	ab179064 - H2A.Z K4ac + K7ac + K11ac
F8	7.0	ab179065 - H2B unmod(1-30)
F9	4.4	ab179066 - H2B K5ac
F10	7.0	ab179067 - H2B K5cr
F11	6.6	ab179068 - H2B K5fo
F12	6.5	ab179069 - H2B P10hyd
F13	5.9	ab179070 - H2B K11me1
F14	4.2	ab179071 - H2B K11me2
F15	6.4	ab179072 - H2B K11me3
F16	5.9	ab179073 - H2B K11cr
F17	4.4	ab179074 - H2B K12ac
F18	3.9	ab179075 - H2B K12cr
F19	4.0	ab179076 - H2B K12me1
F20	11.2	ab179077 - H2B K12me2
F21	1.5	ab179078 - H2B K12me3
F22	32.9	ab179079 - H2B K12ac + K15ac
F23	2.3	ab179080 - H2B S14p
F24	13.8	ab179081 - H2B K15ac
F25	2.9	ab179082 - H2B K15cr
G1	7.7	ab179083 - H2B K15me1
G2	10.2	ab179084 - H2B K15me2
G3	4.8	ab179085 - H2B K15me3
G4	7.6	ab179086 - H2B K16ac
G5	3.4	ab179087 - H2B K16cr
G6	10.1	ab179088 - H2B K20ac
G7	5.1	ab179089 - H2B K20cr
G8	9.4	ab179090 - H2B K20me1
G9	7.5	ab179091 - H2B K20me2
G10	11.4	ab179092 - H2B K20me3
G11	4.5	ab179093 - H2B K23cr
G12	8.3	ab179094 - H2B K23me1
G13	6.8	ab179095 - H2B K23me3
G14	4.2	ab179096 - H2B K24ac
G15	8.0	ab179097 - H2B unmod(26-56)
G16	10.3	ab179098 - H2B K34cr
G17	3.6	ab179099 - H2B K34fo
G18	4.9	ab179100 - H2B Y37hyd
G19	4.5	ab179101 - H2B Y37p
G20	6.5	ab179102 - H2B K43me2
G21	4.2	ab179103 - H2B K43me3
G22	4.2	ab179104 - H2B K46fo
G23	2.9	ab179105 - H2B K46ac
G24	2.6	ab179106 - H2B K46me1
G25	2.7	ab179107 - H2B K46me2
H1	6.8	ab179108 - H2B K46me3
H2	1.7	ab179109 - H2B unmod(51-81)
H3	3.2	ab179110 - H2B K57me1
H4	2.2	ab179111 - H2B K57me2
H5	3.6	ab179112 - H2B K57me3
H6	12.7	ab179113 - H2B unmod(73-103)
H7	5.5	ab179114 - H2B R79me1
H8	9.6	ab179115 - H2B R79me2s
H9	5.7	ab179116 - H2B R79me2a
H10	9.0	ab179117 - H2B Y83hyd
H11	9.2	ab179118 - H2B K85me1
H12	13.5	ab179119 - H2B K85me2
H13	4.8	ab179120 - H2B K85me3
H14	9.0	ab179121 - H2B unmod(94-125)

H15	7.1	ab179122 - H2B R99me1
H16	9.4	ab179123 - H2B R99me2s
H17	9.5	ab179124 - H2B R99me2a
H18	4.1	ab179125 - H2B K108fo
H19	5.6	ab179126 - H2B K116fo
H20	3.7	ab179127 - H2B K116me1
H21	7.0	ab179128 - H2B K116me2
H22	5.4	ab179129 - H2B K116me3
H23	5.9	ab179130 - H2B K116ubiq
H24	3.4	ab179131 - H2B K120ac
H25	5.5	ab179132 - H2B K120fo
I1	9.3	ab179133 - H2B K120ubiq
I2	2.8	ab179134 - H2B K125ubiq
I3	2.6	ab179135 - ScH2B unmod(102-132)
I4	4.4	ab179136 - H3 unmod(1-33)
I5	2.3	ab179137 - H3 R2cit
I6	1611.2	ab179138 - H3 R2cit + K4me2
I7	1.8	ab179139 - H3 R2cit + K4me3
I8	411.4	ab179140 - H3 R2cit + K4me3 + K9ac + K14ac + K18ac
I9	3.1	ab179141 - H3 R2cit + K4ac + K9ac + K14ac + K18ac
I10	2.2	ab179142 - H3 R2me1
I11	1273.8	ab179143 - H3 R2me1 + K4me2
I12	1.5	ab179144 - H3 R2me1 + K4me3
I13	186.4	ab179145 - H3 R2me1 + K4me3 + K9ac + K14ac + K18ac
I14	2.9	ab179146 - H3 R2me2s
I15	1072.6	ab179147 - H3 R2me2s + K4me2
I16	5.2	ab179148 - H3 R2me2s + K4me3
I17	16.1	ab179149 - H3 R2me2s + K4me3 + K9ac + K14ac + K18ac
I18	9.3	ab179150 - H3 R2me2a
I19	1226.3	ab179151 - H3 R2me2a + K4me2
I20	7.1	ab179152 - H3 R2me2a + K4me3
I21	335.4	ab179153 - H3 R2me2a + K4me3 + S10p
I22	6.0	ab179154 - H3 R2me2a + T3p + K4me3
I23	30.1	ab179155 - H3 R2me2a + K4ac + K9ac + K14ac + K18ac
I24	15.0	ab179156 - H3 R2me2a + K4me3 + K9ac + K14ac + K18ac
I25	12.5	ab179157 - H3 R2me2a + T3p + K4me3 + K9ac + K14ac + K18ac
J1	10.6	ab179158 - H3 R2me2a + K4me3 + K9ac + S10p + K14ac + K18ac
J2	9.2	ab179159 - H3 T3p
J3	10.6	ab179160 - H3 T3p + K4me1
J4	321.5	ab179161 - H3 T3p + K4me2
J5	10.0	ab179162 - H3 T3p + K4me3 + K9ac + K14ac + K18ac
J6	2.1	ab179163 - H3 K4ac
J7	74.5	ab179164 - H3 K4ac + K9ac
J8	3.1	ab179165 - H3 K4ac + K14ac
J9	17.6	ab179166 - H3 K4ac + K18ac
J10	2.6	ab179167 - H3 K4ac + K9me3
J11	10.5	ab179168 - H3 K4ac + K9ac + K14ac
J12	3.0	ab179169 - H3 K4ac + K9ac + K18ac
J13	6.7	ab179170 - H3 K4ac + K14ac + K18ac
J14	4.9	ab179171 - H3 K4ac + K9ac + K14ac + K18ac
J15	37.4	ab179172 - H3 K4ac + K9me1 + K14ac + K18ac
J16	86.9	ab179173 - H3 K4ac + K9me2 + K14ac + K18ac
J17	27.6	ab179174 - H3 K4ac + K9me3 + K14ac + K18ac
J18	9.0	ab179175 - H3 K4ac + K9ac + S10p + K14ac + K18ac
J19	31.8	ab179176 - H3 K4cr
J20	297.7	ab179177 - H3 K4me1 + K9me2
J21	32.7	ab179178 - H3 K4me1 + K9ac + K18ac
J22	39.5	ab179179 - H3 K4me1 + K9ac + K14ac + K18ac
J23	60.6	ab179180 - H3 K4me2 + T6p
J24	2188.7	ab179181 - H3 K4me2 + K9me2
J25	4259.1	ab179182 - H3 K4me2 + K9ac + K18ac
K1	2837.6	ab179183 - H3 K4me2 + K9ac + K14ac + K18ac
K2	2.8	ab179184 - H3 K4me3 + T6p
K3	9.8	ab179185 - H3 K4me3 + K9ac

K4	23.3	ab179186 - H3 K4me3 + K9me2
K5	10.2	ab179187 - H3 K4me3 + K9me3
K6	13.1	ab179188 - H3 K4me3 + S10p
K7	5.8	ab179189 - H3 K4me3 + K14ac
K8	8.4	ab179190 - H3 K4me3 + K18ac
K9	4.0	ab179191 - H3 K4me3 + R8me2s + K9me3
K10	6.4	ab179192 - H3 K4me3 + R8me2a + K9me3
K11	101.6	ab179193 - H3 K4me3 + K9ac + S10p
K12	8.8	ab179194 - H3 K4me3 + K9ac + K14ac
K13	35.1	ab179195 - H3 K4me3 + K9ac + K18ac
K14	29.0	ab179196 - H3 K4me3 + K14ac + K18ac
K15	30.7	ab179197 - H3 K4me3 + K9ac + K14ac + K18ac
K16	15.4	ab179198 - H3 K4me3 + K9ac + S10p + K14ac + K18ac
K17	16.1	ab179199 - H3 K4me3 + T6p + K9ac + K14ac + K18ac
K18	916.8	ab179200 - H3 K4me2 + T6p + K9ac + K14ac + K18ac
K19	63.3	ab179201 - H3 K4bio + K9bio + K18bio
K20	20.5	ab179202 - H3 T6p + K9me3
K21	23.9	ab179203 - H3 T6p + R8me2a + K9me3
K22	10.2	ab179204 - H3 T6p + K9ac + K14ac + K18ac
K23	35.3	ab179205 - H3 R8me1
K24	12.8	ab179206 - H3 R8me1 + K9me1
K25	22.1	ab179207 - H3 R8me1 + K9me2
L1	4.0	ab179208 - H3 R8me1 + K9me3
L2	6.2	ab179209 - H3 R8me2s
L3	3.0	ab179210 - H3 R8me2s + K9me1
L4	14.4	ab179211 - H3 R8me2s + K9me2
L5	4.3	ab179212 - H3 R8me2s + K9me3
L6	7.1	ab179213 - H3 R8me2a
L7	7.4	ab179214 - H3 R8me2a + K9me1
L8	16.4	ab179215 - H3 R8me2a + K9me2
L9	7.7	ab179216 - H3 R8me2a + K9me3
L10	7.3	ab179217 - H3 R8cit
L11	8.2	ab179218 - H3 K9ac + K14ac
L12	11.0	ab179219 - H3 K9ac + K18ac
L13	4.4	ab179220 - H3 K9ac + K14ac + K18ac
L14	6.1	ab179221 - H3 K9ac + K14ac + K18ac + K23ac + K27ac
L15	7.6	ab179222 - H3 K9ac + S10p
L16	9.5	ab179223 - H3 K9cr
L17	17.7	ab179224 - H3 K9me2
L18	7.0	ab179225 - H3 K9me2 + S10p
L19	20.1	ab179226 - H3 K9me2 + K27me2
L20	7.0	ab179227 - H3 S10p + S28p
L21	8.1	ab179228 - H3 S10p K14ac
L22	3.9	ab179229 - H3 K14me1
L23	5.3	ab179230 - H3 K14me2
L24	9.5	ab179231 - H3 K14me3
L25	5.8	ab179232 - H3 K14ac + K18ac
M1	5.4	ab179233 - H3 R17cit
M2	6.3	ab179234 - H3 R17me1
M3	5.3	ab179235 - H3 R17me2s
M4	8.6	ab179236 - H3 R17me2a + K18ac
M5	2.4	ab179237 - H3 K18cr
M6	4.4	ab179238 - H3 K18me1
M7	6.3	ab179239 - H3 K18me2
M8	16.2	ab179240 - H3 K18me3 + K36me3
M9	10.1	ab179241 - H3 K23cr
M10	10.5	ab179242 - H3 K23me1
M11	11.5	ab179243 - H3 K23me2
M12	9.7	ab179244 - H3 K23me3
M13	10.5	ab179245 - H3 T22p + K23me3
M14	7.5	ab179246 - H3 K23pr
M15	12.4	ab179247 - H3 unmod(20-51)
M16	6.5	ab179248 - H3 R26me1
M17	6.6	ab179249 - H3 R26me2a

M18	7.0	ab179250 - H3 R26me2a + K27me1
M19	6.8	ab179251 - H3 R26me2a + K27me3
M20	3.5	ab179252 - H3 K27ac + S28p
M21	4.4	ab179253 - H3 K27cr
M22	5.4	ab179254 - H3 K27me1
M23	2.7	ab179255 - H3 K27me2 + S28p
M24	5.8	ab179256 - H3 K27me3
M25	3.7	ab179257 - H3 K27me3 + S28p
N1	6.6	ab179258 - H3 P30hpro
N2	8.9	ab179259 - H3 P30(5,5me2)
N3	6.5	ab179260 - H3 P30
N4	8.9	ab179261 - H3 P30mutA
N5	6.5	ab179262 - H3 T32GlcNAc
N6	6.2	ab179263 - H3 K36ac
N7	10.0	ab179264 - H3 P38hpro
N8	6.8	ab179265 - H3 P38(5,5me2)
N9	12.9	ab179266 - H3 P38
N10	14.5	ab179267 - H3 P38mutA
N11	6.9	ab179268 - H3 Y41p
N12	6.8	ab179269 - H3 unmod(50-80)
N13	5.8	ab179270 - H3 K56ac
N14	4.6	ab179271 - H3 K56cr
N15	4.1	ab179272 - H3 K56fo
N16	7.2	ab179273 - H3 K56me1
N17	8.6	ab179274 - H3 K56me2
N18	5.3	ab179275 - H3 K56me3
N19	9.1	ab179276 - H3 R63me1
N20	6.1	ab179277 - H3 R63me2s
N21	8.8	ab179278 - H3 R63me2a
N22	6.8	ab179279 - H3 K64ac
N23	6.1	ab179280 - H3 K64me2
N24	8.8	ab179281 - H3 K64me3
N25	7.5	ab179282 - H3 unmod(73-103)
O1	7.3	ab179283 - H3 K79ac
O2	4.8	ab179284 - H3 K79fo
O3	3.9	ab179285 - H3 T80p
O4	4.6	ab179286 - H3 K79me3 + T80p
O5	4.6	ab179287 - H3 unmod(99-135)
O6	2.9	ab179288 - H3 T118O-ac
O7	4.4	ab179289 - H3 K122fo
O8	4.0	ab179290 - H3 K122me1
O9	4.4	ab179291 - H3 R128me1
O10	9.2	ab179292 - H3 R128me2s
O11	11.4	ab179293 - H3.3 unmod(10-40)
O12	3.7	ab179294 - H3.3 S28p
O13	8.3	ab179295 - H3.3 S28p + S31p
O14	10.4	ab179296 - H3.3 S31p
O15	11.9	ab179297 - H4 unmod(1-30)
O16	10.2	ab179298 - H4 S1p + R3me1
O17	7.4	ab179299 - H4 S1p + R3me2s
O18	11.1	ab179300 - H4 S1p + R3me2a
O19	9.3	ab179301 - H4 S1p + K5ac + K8ac + K12ac + K16ac
O20	11.1	ab179302 - H4 R3cit
O21	18.9	ab179303 - H4 R3me1 + K5ac + K8ac + K12ac + K16ac + K20ac
O22	8.4	ab179304 - H4 R3me2s
O23	14.2	ab179305 - H4 R3me2s + K5ac + K8ac + K12ac + K16ac + K20ac
O24	10.2	ab179306 - H4 R3me2a
O25	4.5	ab179307 - H4 R3me2a + K5ac
P1	30.3	ab179308 - H4 R3me2a + K5ac + K8ac + K12ac + K16ac + K20ac
P2	3.8	ab179309 - H4 K5ac
P3	11.8	ab179310 - H4 K5ac + K8ac
P4	9.9	ab179311 - H4 K5ac + K12ac
P5	9.0	ab179312 - H4 K5ac + K16ac
P6	7.8	ab179313 - H4 K5ac + K8ac + K12ac

P7	10.3	ab179314 - H4 K5ac + K8me1 + K12ac
P8	6.8	ab179315 - H4 K5ac + K8ac + K12ac + K16ac
P9	8.0	ab179316 - H4 K5ac + K8me1 + K12ac + K16ac
P10	11.7	ab179317 - H4 K5ac + K8ac + K12me1 + K16ac
P11	8.8	ab179318 - H4 K5me1
P12	4.2	ab179319 - H4 K5me1 + K8ac + K12me1
P13	10.8	ab179320 - H4 K5me1 + K8me1 + K12me1
P14	3.6	ab179321 - H4 K5me1 + K8ac + K12ac + K16ac
P15	7.7	ab179322 - H4 K5me2
P16	4.0	ab179323 - H4 K5me3
P17	8.6	ab179324 - H4 K5cr
P18	12.1	ab179325 - H4 K8ac + K12ac
P19	7.0	ab179326 - H4 K8ac + K16ac
P20	9.4	ab179327 - H4 K8cr
P21	9.6	ab179328 - H4 K8me1
P22	9.3	ab179329 - H4 K12ac
P23	8.2	ab179330 - H4 K12ac + K16ac
P24	9.7	ab179331 - H4 K12ac + K16ac + K20me2
P25	6.1	ab179332 - H4 K12ac + K16ac + K20me3
Q1	9.4	ab179333 - H4 K5mutQ + K8mutQ + K12ac + K16mutQ + K20mutQ
Q2	4.4	ab179334 - H4 K12cr
Q3	6.3	ab179335 - H4 K12fo
Q4	6.6	ab179336 - H4 K12me1
Q5	9.4	ab179337 - H4 K12me2
Q6	4.9	ab179338 - H4 K12me3
Q7	6.6	ab179339 - H4 K16ac
Q8	10.1	ab179340 - H4 K16ac + K20ac
Q9	7.7	ab179341 - H4 K16cr
Q10	10.4	ab179342 - H4 K16me1
Q11	7.6	ab179343 - H4 K16me3
Q12	5.1	ab179344 - H4 R17me1
Q13	8.1	ab179345 - H4 R17me2s
Q14	8.9	ab179346 - H4 R17me2a
Q15	6.7	ab179347 - H4 R19me1
Q16	6.8	ab179348 - H4 R19me2s
Q17	11.3	ab179349 - H4 R19me2a
Q18	5.8	ab179350 - H4 K20ac
Q19	7.5	ab179351 - H4 R23me1
Q20	8.0	ab179352 - H4 R23me2s
Q21	6.8	ab179353 - H4 R23me2a
Q22	10.7	ab179354 - H4 unmod(24-54)
Q23	9.3	ab179355 - H4 K31fo
Q24	7.6	ab179356 - H4 K31me1
Q25	7.6	ab179357 - H4 K31me2
R1	5.8	ab179358 - H4 K31me3
R2	10.4	ab179359 - H4 P32hyd
R3	4.3	ab179360 - H4 R35me1
R4	7.7	ab179361 - H4 R35me2s
R5	4.9	ab179362 - H4 R35me2a
R6	6.9	ab179363 - H4 S47O-ac
R7	6.1	ab179364 - H4 S47GlcNAc
R8	9.8	ab179365 - H4 S47p
R9	7.2	ab179366 - H4 unmod(45-75)
R10	7.0	ab179367 - H4 Y51hyd
R11	9.9	ab179368 - H4 R55me1
R12	9.4	ab179369 - H4 R55me2s
R13	7.9	ab179370 - H4 R55me2a
R14	8.5	ab179371 - H4 K59fo
R15	8.7	ab179372 - H4 K59me1
R16	5.8	ab179373 - H4 K59me2
R17	4.2	ab179374 - H4 R67me1
R18	4.1	ab179375 - H4 R67me2s
R19	4.6	ab179376 - H4 R67me2a
R20	7.4	ab179377 - H4 unmod(71-102)

R21	4.6	ab179378 - H4 K77me1
R22	4.9	ab179379 - H4 K77me2
R23	8.8	ab179380 - H4 K77me3
R24	9.8	ab179381 - H4 K79fo
R25	8.2	ab179382 - H4 K79me1
S1	9.1	ab179383 - H4 K79me3
S2	7.4	ab179384 - H4 T82O-ac
S3	10.4	ab179385 - H4 Y88hyd
S4	6.2	ab179386 - H4 Y88nitro
S5	7.1	ab179387 - H4 K91fo
S6	10.2	ab179403 - H2A K5ac
S7	9.6	ab179404 - H2A K9me1
S8	5.7	ab18504 - H2B S32p
S9	5.7	ab19828 - ScH2A S129p
S10	5.6	ab20631 - H3 R26cit
S11	7.3	ab20635 - H3 T3p + K4me3
S12	6.8	ab21043 - H3 K37me1
S13	1.9	ab21998 - H1.4 K25me2
S14	4.2	ab222079 - H4 K59me3
S15	4.4	ab222080 - H3 K122me2
S16	8.0	ab222082 - H3 K122me3
S17	25.2	ab22214 - H3 K37me2
S18	6.5	ab22398 - H2A R3me2a
S19	7.3	ab22399 - H2A R3me2s
S20	9.3	ab23539 - H3 T22p
S21	4.2	ab24003 - H3 K18ac
S22	6.1	ab24404 - H3 K27ac
S23	6.6	ab24417 - H3 K37me3
S24	8.9	ab24444 - H3 T11p
S25	6.0	ab24659 - H2A K127me3
T1	5.7	ab24660 - H2A K125me3
T2	7.7	ab26307 - H3 T45p
T3	5.0	ab26342 - H2B K5me2
T4	8.6	ab26739 - H4 T30p
T5	6.8	ab27896 - H4 T96p
T6	10.0	ab27898 - H4 T80p
T7	5.7	ab28756 - H2A R3me1
T8	6.8	ab32876 - H3 R2cit + R8cit + R17cit
T9	6.9	ab33047 - H2A K127me2
T10	3.6	ab33505 - H3 T118p
T11	4.8	ab33594 - H3 R128me2a
T12	3.7	ab34466 - H3 K122ac
T13	3.8	ab36927 - H2A S122p
T14	8.5	ab4555 - H3 K79me1
T15	6.0	ab4556 - H3 K79me2
T16	6.0	ab4557 - H3 K79me3
T17	12.3	ab4560 - H4 K79me2
T18	10.0	ab46854 - H3 K18me3
T19	11.7	ab48359 - H3 K23ac
T20	12.1	ab53244 - H3 K64me1
T21	8.9	ab54016 - H2A K9ac
T22	8.8	ab5499 - H3 S28p
T23	1474.7	ab7768 - H3 K4me2
T24	4.0	ab89018 - H1.4 S26p
T25	817.8	ab92374 - H3 K4me3 (HEX1)
U1	9.3	ab95444 - H4 K16me2