


Anti-Histone H3 (tri methyl K9) antibody [EPR16601] – ChIP Grade (ab176916)


Application	ChIP-seq
Target	Histone H3 (tri methyl K9)
Sample	MEF
Treatment	None
Antibody amount	4 µg
Cell number	1X10 ⁷
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	mm39
Reads	30 million

Chromosome 17




Copyright (c) 2023 Abcam plc.

Chromosome 7



Copyright (c) 2023 Abcam plc.

Chromosome 7



Chromosome 7

