






Anti-Histone H3 (mono+di+tri methyl K79) antibody [EPR17468] - ChIP Grade (ab177185)

Target	Histone H3 (mono+di+tri methyl K79)
Sample	HeLa
Treatment	None
Antibody amount	4 µg
Cell number	1×10 ⁷
Sequencing	NovaSeq 6000 (PE 150)
Peak calling	MACS 2.1.0
Reference genome	hg38






Copyright (c) 2021 Abcam plc.



Copyright (c) 2021 Abcam plc.



Copyright (c) 2021 Abcam plc.