

Peptide array: H3 R26me2s (ab177284)

Area under curve (AUC) data for all peptides on the pad. Antigen-containing peptides are highlighted in green.

Cell	AUC	Peptide
A1	5.1	ab10112 - H4 K91ac
A2	5.5	ab10121 - H2B K85ac
A3	4.4	ab10127 - H2A K95me2
A4	4.5	ab10138 - H1.4 T17p
A5	4.4	ab10139 - H1.4 T146p
A6	6.0	ab113997 - H2A.X Y142p
A7	3.4	ab11477 - H3 S10p
A8	6.4	ab12946 - H2B K5me3
A9	3.4	ab12949 - H2B K23me2
A10	4.9	ab12952 - H2A K127me1
A11	5.9	ab13211 - H2B K5me1
A12	3.9	ab1340 - H3 K4me1
A13	4.4	ab13570 - H2B K43me1
A14	3.4	ab13835 - HA
A15	6.3	ab13837 - Myc
A16	5.4	ab140408 - ScH2B T129p
A17	3.6	ab14103 - H3 T6p
A18	4.0	ab14724 - H4 S1p
A19	5.3	ab14799 - H3 T32p
A20	4.4	ab14943 - His
A21	3.3	ab14952 - H2A R77me2s
A22	5.3	ab14964 - H4 K20me2
A23	4.9	ab14997 - ScH2A K21ac
A24	9283.3	ab154782 - H3 R26me2s
A25	4.7	ab15591 - H3 K14ac
B1	2.5	ab15644 - H3 K9me3 + S10p
B2	4.0	ab15645 - H2A.X S139p
B3	2.9	ab15824 - H4 K8ac
B4	3.3	ab15829 - V5
B5	2.7	ab16635 - H3 K9ac
B6	7.5	ab166685 - H2B S112GlcNAc
B7	5.4	ab166688 - H2B S36GlcNAc
B8	4.3	ab16935 - H3 R17me2a
B9	4.7	ab17043 - H4 K20me1
B10	4.2	ab17162 - H3 P16hyd + P30hyd
B11	3.9	ab17567 - H4 K20me3
B12	3.2	ab17587 - H1.4 K25me3
B13	5.6	ab175905 - H2B K11ac
B14	3.8	ab176211 - H2B K23ac
B15	3.1	ab17632 - H4 K5bio
B16	4.7	ab1771 - H3 K9me1
B17	3.6	ab1773 - H3 K9me3
B18	3.7	ab17770 - H4 R3me1
B19	6.2	ab1781 - H3 K27me2
B20	99.5	ab1783 - H3 K36me1
B21	4.1	ab1784 - H3 K36me2
B22	5.3	ab1785 - H3 K36me3
B23	3.6	ab178667 - ScH2A unmod(1-30)
B24	4.7	ab178668 - ScH2A unmod(102-132)
B25	5.1	ab178670 - H2A.Z unmod(1-30)
C1	7.8	ab178833 - Flag
C2	3.0	ab178835 - H2B R99mutO
C3	5.8	ab178836 - H4 R23mutO
C4	3.8	ab178837 - H4 G48mutD
C5	12.3	ab178987 - H1.4 unmod(1-30)
C6	2.9	ab178988 - H1.4 S1p
C7	6.4	ab178989 - H1.4 T3p
C8	5.2	ab178990 - H1.4 S1p + T3p
C9	2.2	ab178991 - H1.4 K25me1
C10	3.3	ab178992 - H1.4 K25me3 + S26p
C11	3.4	ab178993 - H1.4 unmod(130-160)

C12	4.7	ab178994 - H2A unmod(1-30)
C13	3.2	ab178995 - H2A S1p
C14	4.3	ab178996 - H2A S1p + K5ac
C15	3.0	ab178997 - H2A S1p + R3me2a + K5ac
C16	3.9	ab178998 - H2A S1p + R3me2s + K5ac
C17	3.7	ab178999 - H2A S1p + R3cit + K5ac
C18	3.1	ab179000 - H2A S1ac + S1p + K5ac + K9ac + K13ac + K15ac
C19	3.8	ab179001 - H2A S1p + K5ac + K9ac + K13ac + K15ac
C20	2.4	ab179002 - H2A R3cit
C21	4.5	ab179003 - H2A R3cit + K5ac
C22	4.3	ab179004 - H2A R3me2a + K5ac
C23	9.3	ab179005 - H2A R3me2s + K5ac
C24	4.4	ab179006 - H2A S1ac + K5ac + K9ac + K13ac + K15ac
C25	6.1	ab179007 - H2A K5ac + K9ac + K13ac + K15ac
D1	4.8	ab179008 - H2A K9me2
D2	2.5	ab179009 - H2A K9me3
D3	3.9	ab179010 - H2A R11me1
D4	2.9	ab179011 - H2A R11me2s
D5	5.1	ab179012 - H2A R11me2a
D6	3.6	ab179013 - H2A K13ac
D7	4.6	ab179014 - H2A K15ac
D8	3.6	ab179015 - H2A R17me1
D9	9.9	ab179016 - H2A R17me2s
D10	4.4	ab179017 - H2A R17me2a
D11	5.5	ab179018 - H2A unmod(20-50)
D12	4.0	ab179019 - H2A P26hyd
D13	2.5	ab179020 - H2A R29me1
D14	8.9	ab179021 - H2A R29me2s
D15	5.7	ab179022 - H2A R29me2a
D16	3.3	ab179023 - H2A K36ac
D17	3.2	ab179024 - H2A K36cr
D18	5.2	ab179025 - H2A Y39hyd
D19	4.0	ab179026 - H2A R42me1
D20	54.7	ab179027 - H2A R42me2s
D21	2.7	ab179028 - H2A R42me2a
D22	2.3	ab179029 - H2A unmod(71-103)
D23	4.2	ab179030 - H2A R77me1
D24	6.8	ab179031 - H2A R77me2a
D25	3.5	ab179032 - H2A R88me1
E1	4.5	ab179033 - H2A R88me2s
E2	4.1	ab179034 - H2A R88me2a
E3	4.8	ab179035 - H2A K95me1
E4	4.5	ab179036 - H2A K95me3
E5	3.3	ab179037 - H2A unmod(92-123)
E6	3.1	ab179038 - H2A T101GlcNAc
E7	4.1	ab179039 - H2A Q104me1
E8	49.2	ab179040 - H2A unmod(112-129)
E9	3.5	ab179041 - H2A K118me1
E10	3.5	ab179042 - H2A K118me2
E11	3.4	ab179043 - H2A K118me3
E12	3.6	ab179044 - H2A K118fo
E13	3.3	ab179045 - H2A K118cr
E14	3.6	ab179046 - H2A K119cr
E15	3.3	ab179047 - H2A K118ubiq
E16	4.4	ab179048 - H2A K119ubiq
E17	3.4	ab179049 - H2A T120p
E18	3.3	ab179050 - H2A K125me1
E19	3.9	ab179051 - H2A K125me2
E20	2.7	ab179052 - H2A K125cr
E21	8.0	ab179053 - H2A.X unmod(1-30)
E22	4.7	ab179054 - H2A.X K5ac
E23	3.6	ab179055 - H2A.X unmod(114-142)
E24	4.8	ab179056 - H2A.X K118ubiq
E25	6.0	ab179057 - H2A.X K119ubiq

F1	4.7	ab179058 - H2A.X T120p
F2	5.8	ab179059 - H2A.X K127ubiq
F3	4.1	ab179060 - H2A.X S139p + Y142p
F4	4.9	ab179061 - H2A.Z K4ac
F5	4.0	ab179062 - H2A.Z K7ac
F6	3.0	ab179063 - H2A.Z K11ac
F7	3.2	ab179064 - H2A.Z K4ac + K7ac + K11ac
F8	5.2	ab179065 - H2B unmod(1-30)
F9	3.3	ab179066 - H2B K5ac
F10	6.0	ab179067 - H2B K5cr
F11	4.9	ab179068 - H2B K5fo
F12	4.8	ab179069 - H2B P10hyd
F13	15.0	ab179070 - H2B K11me1
F14	3.3	ab179071 - H2B K11me2
F15	4.7	ab179072 - H2B K11me3
F16	4.0	ab179073 - H2B K11cr
F17	4.7	ab179074 - H2B K12ac
F18	3.5	ab179075 - H2B K12cr
F19	3.0	ab179076 - H2B K12me1
F20	3.7	ab179077 - H2B K12me2
F21	4.6	ab179078 - H2B K12me3
F22	7.3	ab179079 - H2B K12ac + K15ac
F23	5.7	ab179080 - H2B S14p
F24	4.5	ab179081 - H2B K15ac
F25	4.9	ab179082 - H2B K15cr
G1	3.9	ab179083 - H2B K15me1
G2	4.1	ab179084 - H2B K15me2
G3	4.4	ab179085 - H2B K15me3
G4	4.1	ab179086 - H2B K16ac
G5	5.7	ab179087 - H2B K16cr
G6	5.5	ab179088 - H2B K20ac
G7	1.4	ab179089 - H2B K20cr
G8	5.4	ab179090 - H2B K20me1
G9	4.0	ab179091 - H2B K20me2
G10	4.0	ab179092 - H2B K20me3
G11	3.7	ab179093 - H2B K23cr
G12	6.8	ab179094 - H2B K23me1
G13	4.4	ab179095 - H2B K23me3
G14	4.4	ab179096 - H2B K24ac
G15	7.6	ab179097 - H2B unmod(26-56)
G16	6.4	ab179098 - H2B K34cr
G17	4.4	ab179099 - H2B K34fo
G18	4.8	ab179100 - H2B Y37hyd
G19	4.5	ab179101 - H2B Y37p
G20	5.4	ab179102 - H2B K43me2
G21	4.7	ab179103 - H2B K43me3
G22	3.3	ab179104 - H2B K46fo
G23	4.1	ab179105 - H2B K46ac
G24	4.6	ab179106 - H2B K46me1
G25	4.5	ab179107 - H2B K46me2
H1	4.3	ab179108 - H2B K46me3
H2	4.6	ab179109 - H2B unmod(51-81)
H3	4.7	ab179110 - H2B K57me1
H4	3.7	ab179111 - H2B K57me2
H5	4.7	ab179112 - H2B K57me3
H6	4.5	ab179113 - H2B unmod(73-103)
H7	4.4	ab179114 - H2B R79me1
H8	5.9	ab179115 - H2B R79me2s
H9	2.8	ab179116 - H2B R79me2a
H10	3.9	ab179117 - H2B Y83hyd
H11	4.2	ab179118 - H2B K85me1
H12	5.2	ab179119 - H2B K85me2
H13	5.0	ab179120 - H2B K85me3
H14	5.7	ab179121 - H2B unmod(94-125)

H15	3.3	ab179122 - H2B R99me1
H16	6.2	ab179123 - H2B R99me2s
H17	4.5	ab179124 - H2B R99me2a
H18	4.0	ab179125 - H2B K108fo
H19	3.4	ab179126 - H2B K116fo
H20	4.5	ab179127 - H2B K116me1
H21	3.6	ab179128 - H2B K116me2
H22	3.5	ab179129 - H2B K116me3
H23	4.0	ab179130 - H2B K116ubiq
H24	4.0	ab179131 - H2B K120ac
H25	4.7	ab179132 - H2B K120fo
I1	3.4	ab179133 - H2B K120ubiq
I2	5.2	ab179134 - H2B K125ubiq
I3	4.2	ab179135 - ScH2B unmod(102-132)
I4	44.5	ab179136 - H3 unmod(1-33)
I5	3.1	ab179137 - H3 R2cit
I6	5.5	ab179138 - H3 R2cit + K4me2
I7	3.3	ab179139 - H3 R2cit + K4me3
I8	7.3	ab179140 - H3 R2cit + K4me3 + K9ac + K14ac + K18ac
I9	4.7	ab179141 - H3 R2cit + K4ac + K9ac + K14ac + K18ac
I10	3.5	ab179142 - H3 R2me1
I11	4.0	ab179143 - H3 R2me1 + K4me2
I12	3.9	ab179144 - H3 R2me1 + K4me3
I13	2.1	ab179145 - H3 R2me1 + K4me3 + K9ac + K14ac + K18ac
I14	5.6	ab179146 - H3 R2me2s
I15	8.9	ab179147 - H3 R2me2s + K4me2
I16	34.9	ab179148 - H3 R2me2s + K4me3
I17	24.2	ab179149 - H3 R2me2s + K4me3 + K9ac + K14ac + K18ac
I18	3.0	ab179150 - H3 R2me2a
I19	4.5	ab179151 - H3 R2me2a + K4me2
I20	3.4	ab179152 - H3 R2me2a + K4me3
I21	2.2	ab179153 - H3 R2me2a + K4me3 + S10p
I22	3.0	ab179154 - H3 R2me2a + T3p + K4me3
I23	3.5	ab179155 - H3 R2me2a + K4ac + K9ac + K14ac + K18ac
I24	4.1	ab179156 - H3 R2me2a + K4me3 + K9ac + K14ac + K18ac
I25	4.5	ab179157 - H3 R2me2a + T3p + K4me3 + K9ac + K14ac + K18ac
J1	4.1	ab179158 - H3 R2me2a + K4me3 + K9ac + S10p + K14ac + K18ac
J2	3.2	ab179159 - H3 T3p
J3	5.5	ab179160 - H3 T3p + K4me1
J4	2.7	ab179161 - H3 T3p + K4me2
J5	3.5	ab179162 - H3 T3p + K4me3 + K9ac + K14ac + K18ac
J6	4.3	ab179163 - H3 K4ac
J7	3.3	ab179164 - H3 K4ac + K9ac
J8	5.6	ab179165 - H3 K4ac + K14ac
J9	3.5	ab179166 - H3 K4ac + K18ac
J10	5.2	ab179167 - H3 K4ac + K9me3
J11	3.8	ab179168 - H3 K4ac + K9ac + K14ac
J12	4.1	ab179169 - H3 K4ac + K9ac + K18ac
J13	3.9	ab179170 - H3 K4ac + K14ac + K18ac
J14	4.1	ab179171 - H3 K4ac + K9ac + K14ac + K18ac
J15	5.3	ab179172 - H3 K4ac + K9me1 + K14ac + K18ac
J16	3.3	ab179173 - H3 K4ac + K9me2 + K14ac + K18ac
J17	2.5	ab179174 - H3 K4ac + K9me3 + K14ac + K18ac
J18	4.5	ab179175 - H3 K4ac + K9ac + S10p + K14ac + K18ac
J19	3.1	ab179176 - H3 K4cr
J20	3.6	ab179177 - H3 K4me1 + K9me2
J21	2.8	ab179178 - H3 K4me1 + K9ac + K18ac
J22	4.8	ab179179 - H3 K4me1 + K9ac + K14ac + K18ac
J23	5.2	ab179180 - H3 K4me2 + T6p
J24	4.8	ab179181 - H3 K4me2 + K9me2
J25	5.1	ab179182 - H3 K4me2 + K9ac + K18ac
K1	2.8	ab179183 - H3 K4me2 + K9ac + K14ac + K18ac
K2	5.1	ab179184 - H3 K4me3 + T6p
K3	2.3	ab179185 - H3 K4me3 + K9ac

K4	4.9	ab179186 - H3 K4me3 + K9me2
K5	4.0	ab179187 - H3 K4me3 + K9me3
K6	2.6	ab179188 - H3 K4me3 + S10p
K7	3.0	ab179189 - H3 K4me3 + K14ac
K8	5.8	ab179190 - H3 K4me3 + K18ac
K9	20.2	ab179191 - H3 K4me3 + R8me2s + K9me3
K10	3.1	ab179192 - H3 K4me3 + R8me2a + K9me3
K11	6.9	ab179193 - H3 K4me3 + K9ac + S10p
K12	2.8	ab179194 - H3 K4me3 + K9ac + K14ac
K13	4.0	ab179195 - H3 K4me3 + K9ac + K18ac
K14	6.8	ab179196 - H3 K4me3 + K14ac + K18ac
K15	2.7	ab179197 - H3 K4me3 + K9ac + K14ac + K18ac
K16	4.4	ab179198 - H3 K4me3 + K9ac + S10p + K14ac + K18ac
K17	4.6	ab179199 - H3 K4me3 + T6p + K9ac + K14ac + K18ac
K18	3.9	ab179200 - H3 K4me2 + T6p + K9ac + K14ac + K18ac
K19	6.7	ab179201 - H3 K4bio + K9bio + K18bio
K20	5.3	ab179202 - H3 T6p + K9me3
K21	15.9	ab179203 - H3 T6p + R8me2a + K9me3
K22	3.9	ab179204 - H3 T6p + K9ac + K14ac + K18ac
K23	4.0	ab179205 - H3 R8me1
K24	2.8	ab179206 - H3 R8me1 + K9me1
K25	6.8	ab179207 - H3 R8me1 + K9me2
L1	3.8	ab179208 - H3 R8me1 + K9me3
L2	78.5	ab179209 - H3 R8me2s
L3	59.3	ab179210 - H3 R8me2s + K9me1
L4	134.6	ab179211 - H3 R8me2s + K9me2
L5	34.5	ab179212 - H3 R8me2s + K9me3
L6	3.3	ab179213 - H3 R8me2a
L7	5.5	ab179214 - H3 R8me2a + K9me1
L8	4.4	ab179215 - H3 R8me2a + K9me2
L9	4.2	ab179216 - H3 R8me2a + K9me3
L10	2.7	ab179217 - H3 R8cit
L11	4.2	ab179218 - H3 K9ac + K14ac
L12	6.9	ab179219 - H3 K9ac + K18ac
L13	3.6	ab179220 - H3 K9ac + K14ac + K18ac
L14	20.1	ab179221 - H3 K9ac + K14ac + K18ac + K23ac + K27ac
L15	4.2	ab179222 - H3 K9ac + S10p
L16	4.8	ab179223 - H3 K9cr
L17	5.6	ab179224 - H3 K9me2
L18	4.9	ab179225 - H3 K9me2 + S10p
L19	6.1	ab179226 - H3 K9me2 + K27me2
L20	4.7	ab179227 - H3 S10p + S28p
L21	5.3	ab179228 - H3 S10p K14ac
L22	4.4	ab179229 - H3 K14me1
L23	3.7	ab179230 - H3 K14me2
L24	5.3	ab179231 - H3 K14me3
L25	4.6	ab179232 - H3 K14ac + K18ac
M1	4.3	ab179233 - H3 R17cit
M2	4.0	ab179234 - H3 R17me1
M3	4.5	ab179235 - H3 R17me2s
M4	4.8	ab179236 - H3 R17me2a + K18ac
M5	3.3	ab179237 - H3 K18cr
M6	6.0	ab179238 - H3 K18me1
M7	4.4	ab179239 - H3 K18me2
M8	51.3	ab179240 - H3 K18me3 + K36me3
M9	4.2	ab179241 - H3 K23cr
M10	4.8	ab179242 - H3 K23me1
M11	3.0	ab179243 - H3 K23me2
M12	5.0	ab179244 - H3 K23me3
M13	4.7	ab179245 - H3 T22p + K23me3
M14	3.4	ab179246 - H3 K23pr
M15	13.2	ab179247 - H3 unmod(20-51)
M16	938.6	ab179248 - H3 R26me1
M17	60.9	ab179249 - H3 R26me2a

M18	1666.5	ab179250 - H3 R26me2a + K27me1
M19	5.8	ab179251 - H3 R26me2a + K27me3
M20	5.2	ab179252 - H3 K27ac + S28p
M21	3.9	ab179253 - H3 K27cr
M22	4.8	ab179254 - H3 K27me1
M23	7.0	ab179255 - H3 K27me2 + S28p
M24	6.0	ab179256 - H3 K27me3
M25	4.7	ab179257 - H3 K27me3 + S28p
N1	4.6	ab179258 - H3 P30hpro
N2	5.7	ab179259 - H3 P30(5,5me2)
N3	3.8	ab179260 - H3 P30
N4	2.4	ab179261 - H3 P30mutA
N5	3.9	ab179262 - H3 T32GlcNAc
N6	3.8	ab179263 - H3 K36ac
N7	5.2	ab179264 - H3 P38hpro
N8	5.0	ab179265 - H3 P38(5,5me2)
N9	5.3	ab179266 - H3 P38
N10	5.3	ab179267 - H3 P38mutA
N11	4.3	ab179268 - H3 Y41p
N12	7.2	ab179269 - H3 unmod(50-80)
N13	3.9	ab179270 - H3 K56ac
N14	3.5	ab179271 - H3 K56cr
N15	4.0	ab179272 - H3 K56fo
N16	3.8	ab179273 - H3 K56me1
N17	6.4	ab179274 - H3 K56me2
N18	4.6	ab179275 - H3 K56me3
N19	3.6	ab179276 - H3 R63me1
N20	4.7	ab179277 - H3 R63me2s
N21	4.0	ab179278 - H3 R63me2a
N22	4.8	ab179279 - H3 K64ac
N23	4.6	ab179280 - H3 K64me2
N24	4.8	ab179281 - H3 K64me3
N25	3.5	ab179282 - H3 unmod(73-103)
O1	5.5	ab179283 - H3 K79ac
O2	4.0	ab179284 - H3 K79fo
O3	3.3	ab179285 - H3 T80p
O4	8.4	ab179286 - H3 K79me3 + T80p
O5	4.1	ab179287 - H3 unmod(99-135)
O6	5.5	ab179288 - H3 T118O-ac
O7	2.4	ab179289 - H3 K122fo
O8	4.7	ab179290 - H3 K122me1
O9	5.2	ab179291 - H3 R128me1
O10	3.5	ab179292 - H3 R128me2s
O11	4.8	ab179293 - H3.3 unmod(10-40)
O12	5.7	ab179294 - H3.3 S28p
O13	4.5	ab179295 - H3.3 S28p + S31p
O14	3.3	ab179296 - H3.3 S31p
O15	8.8	ab179297 - H4 unmod(1-30)
O16	6.6	ab179298 - H4 S1p + R3me1
O17	5.9	ab179299 - H4 S1p + R3me2s
O18	3.5	ab179300 - H4 S1p + R3me2a
O19	3.5	ab179301 - H4 S1p + K5ac + K8ac + K12ac + K16ac
O20	6.5	ab179302 - H4 R3cit
O21	5.8	ab179303 - H4 R3me1 + K5ac + K8ac + K12ac + K16ac + K20ac
O22	4.8	ab179304 - H4 R3me2s
O23	7.6	ab179305 - H4 R3me2s + K5ac + K8ac + K12ac + K16ac + K20ac
O24	2.9	ab179306 - H4 R3me2a
O25	7.6	ab179307 - H4 R3me2a + K5ac
P1	14.1	ab179308 - H4 R3me2a + K5ac + K8ac + K12ac + K16ac + K20ac
P2	4.7	ab179309 - H4 K5ac
P3	7.1	ab179310 - H4 K5ac + K8ac
P4	3.7	ab179311 - H4 K5ac + K12ac
P5	4.2	ab179312 - H4 K5ac + K16ac
P6	3.8	ab179313 - H4 K5ac + K8ac + K12ac

P7	9.3	ab179314 - H4 K5ac + K8me1 + K12ac
P8	5.7	ab179315 - H4 K5ac + K8ac + K12ac + K16ac
P9	14.3	ab179316 - H4 K5ac + K8me1 + K12ac + K16ac
P10	5.3	ab179317 - H4 K5ac + K8ac + K12me1 + K16ac
P11	4.2	ab179318 - H4 K5me1
P12	4.2	ab179319 - H4 K5me1 + K8ac + K12me1
P13	3.6	ab179320 - H4 K5me1 + K8me1 + K12me1
P14	6.5	ab179321 - H4 K5me1 + K8ac + K12ac + K16ac
P15	3.5	ab179322 - H4 K5me2
P16	4.3	ab179323 - H4 K5me3
P17	2.7	ab179324 - H4 K5cr
P18	6.2	ab179325 - H4 K8ac + K12ac
P19	5.6	ab179326 - H4 K8ac + K16ac
P20	7.0	ab179327 - H4 K8cr
P21	4.8	ab179328 - H4 K8me1
P22	4.2	ab179329 - H4 K12ac
P23	3.3	ab179330 - H4 K12ac + K16ac
P24	7.7	ab179331 - H4 K12ac + K16ac + K20me2
P25	4.7	ab179332 - H4 K12ac + K16ac + K20me3
Q1	4.5	ab179333 - H4 K5mutQ + K8mutQ + K12ac + K16mutQ + K20mutQ
Q2	3.8	ab179334 - H4 K12cr
Q3	2.9	ab179335 - H4 K12fo
Q4	5.3	ab179336 - H4 K12me1
Q5	2.0	ab179337 - H4 K12me2
Q6	4.2	ab179338 - H4 K12me3
Q7	3.1	ab179339 - H4 K16ac
Q8	4.5	ab179340 - H4 K16ac + K20ac
Q9	6.0	ab179341 - H4 K16cr
Q10	5.1	ab179342 - H4 K16me1
Q11	5.1	ab179343 - H4 K16me3
Q12	8.0	ab179344 - H4 R17me1
Q13	99.7	ab179345 - H4 R17me2s
Q14	9.8	ab179346 - H4 R17me2a
Q15	3.7	ab179347 - H4 R19me1
Q16	4.3	ab179348 - H4 R19me2s
Q17	4.8	ab179349 - H4 R19me2a
Q18	4.9	ab179350 - H4 K20ac
Q19	5.0	ab179351 - H4 R23me1
Q20	4.7	ab179352 - H4 R23me2s
Q21	6.3	ab179353 - H4 R23me2a
Q22	4.9	ab179354 - H4 unmod(24-54)
Q23	4.4	ab179355 - H4 K31fo
Q24	5.5	ab179356 - H4 K31me1
Q25	4.0	ab179357 - H4 K31me2
R1	5.3	ab179358 - H4 K31me3
R2	5.5	ab179359 - H4 P32hyd
R3	7.2	ab179360 - H4 R35me1
R4	567.4	ab179361 - H4 R35me2s
R5	5.3	ab179362 - H4 R35me2a
R6	5.2	ab179363 - H4 S47O-ac
R7	4.7	ab179364 - H4 S47GlcNAc
R8	3.6	ab179365 - H4 S47p
R9	5.5	ab179366 - H4 unmod(45-75)
R10	5.6	ab179367 - H4 Y51hyd
R11	5.1	ab179368 - H4 R55me1
R12	3.7	ab179369 - H4 R55me2s
R13	5.3	ab179370 - H4 R55me2a
R14	5.1	ab179371 - H4 K59fo
R15	4.6	ab179372 - H4 K59me1
R16	4.2	ab179373 - H4 K59me2
R17	4.9	ab179374 - H4 R67me1
R18	5.4	ab179375 - H4 R67me2s
R19	4.6	ab179376 - H4 R67me2a
R20	6.6	ab179377 - H4 unmod(71-102)

R21	3.7	ab179378 - H4 K77me1
R22	3.8	ab179379 - H4 K77me2
R23	4.2	ab179380 - H4 K77me3
R24	4.6	ab179381 - H4 K79fo
R25	4.5	ab179382 - H4 K79me1
S1	3.7	ab179383 - H4 K79me3
S2	4.5	ab179384 - H4 T82O-ac
S3	5.2	ab179385 - H4 Y88hyd
S4	7.1	ab179386 - H4 Y88nitro
S5	5.5	ab179387 - H4 K91fo
S6	4.1	ab179403 - H2A K5ac
S7	5.0	ab179404 - H2A K9me1
S8	4.2	ab18504 - H2B S32p
S9	4.6	ab19828 - ScH2A S129p
S10	2.5	ab20631 - H3 R26cit
S11	4.9	ab20635 - H3 T3p + K4me3
S12	4.8	ab21043 - H3 K37me1
S13	5.0	ab21998 - H1.4 K25me2
S14	3.3	ab222079 - H4 K59me3
S15	3.5	ab222080 - H3 K122me2
S16	3.6	ab222082 - H3 K122me3
S17	4.5	ab22214 - H3 K37me2
S18	3.3	ab22398 - H2A R3me2a
S19	2.6	ab22399 - H2A R3me2s
S20	4.4	ab23539 - H3 T22p
S21	4.2	ab24003 - H3 K18ac
S22	22.1	ab24404 - H3 K27ac
S23	3.3	ab24417 - H3 K37me3
S24	3.7	ab24444 - H3 T11p
S25	2.8	ab24659 - H2A K127me3
T1	4.0	ab24660 - H2A K125me3
T2	4.9	ab26307 - H3 T45p
T3	4.6	ab26342 - H2B K5me2
T4	5.1	ab26739 - H4 T30p
T5	6.0	ab27896 - H4 T96p
T6	4.2	ab27898 - H4 T80p
T7	3.5	ab28756 - H2A R3me1
T8	4.4	ab32876 - H3 R2cit + R8cit + R17cit
T9	4.9	ab33047 - H2A K127me2
T10	4.4	ab33505 - H3 T118p
T11	3.9	ab33594 - H3 R128me2a
T12	4.1	ab34466 - H3 K122ac
T13	3.6	ab36927 - H2A S122p
T14	5.0	ab4555 - H3 K79me1
T15	3.8	ab4556 - H3 K79me2
T16	3.7	ab4557 - H3 K79me3
T17	4.4	ab4560 - H4 K79me2
T18	4.1	ab46854 - H3 K18me3
T19	4.8	ab48359 - H3 K23ac
T20	4.9	ab53244 - H3 K64me1
T21	4.4	ab54016 - H2A K9ac
T22	4.9	ab5499 - H3 S28p
T23	2.7	ab7768 - H3 K4me2
T24	2.6	ab89018 - H1.4 S26p
T25	4.3	ab92374 - H3 K4me3 (HEX1)
U1	3.9	ab95444 - H4 K16me2