


Anti-Histone H2A antibody [ERP17470] - ChIP Grade (ab177308)


Target	Histone H2A
Sample	HeLa
Treatment	None
Antibody amount	2ug
Cell number	2x10^5
Sequencing	NovaSeq 6000 (PE 150)
Peak calling	MACS 2.1.0
Reference genome	hg38

Chromosome 7




Copyright (c) 2021 Abcam plc.

Chromosome 12




Copyright (c) 2021 Abcam plc.

Chromosome 1



Copyright (c) 2021 Abcam plc.

Chromosome 11



Copyright (c) 2021 Abcam plc.