

Peptide array: H2A K9Ac (ab177312)

Area under curve (AUC) data for all peptides on the pad. Antigen-containing peptides are highlighted in green.

Cell	AUC	Peptide
A1	6.1	ab10112 - H4 K91ac
A2	4.1	ab10121 - H2B K85ac
A3	34.2	ab10127 - H2A K95me2
A4	4.1	ab10138 - H1.4 T17p
A5	3.9	ab10139 - H1.4 T146p
A6	46.6	ab113997 - H2A.X Y142p
A7	4.2	ab11477 - H3 S10p
A8	3.2	ab12946 - H2B K5me3
A9	3.9	ab12949 - H2B K23me2
A10	2.4	ab12952 - H2A K127me1
A11	4.4	ab13211 - H2B K5me1
A12	9.0	ab1340 - H3 K4me1
A13	5.6	ab13570 - H2B K43me1
A14	7.2	ab13835 - HA
A15	6.2	ab13837 - Myc
A16	12.9	ab140408 - ScH2B T129p
A17	5.6	ab14103 - H3 T6p
A18	4.1	ab14724 - H4 S1p
A19	3.8	ab14799 - H3 T32p
A20	6.3	ab14943 - His
A21	7.2	ab14952 - H2A R77me2s
A22	5.1	ab14964 - H4 K20me2
A23	5.3	ab14997 - ScH2A K21ac
A24	3.9	ab154782 - H3 R26me2s
A25	3.3	ab15591 - H3 K14ac
B1	3.5	ab15644 - H3 K9me3 + S10p
B2	22.3	ab15645 - H2A.X S139p
B3	7.1	ab15824 - H4 K8ac
B4	5.7	ab15829 - V5
B5	12.3	ab16635 - H3 K9ac
B6	4.0	ab166685 - H2B S112GlcNAc
B7	4.4	ab166688 - H2B S36GlcNAc
B8	4.4	ab16935 - H3 R17me2a
B9	7.1	ab17043 - H4 K20me1
B10	3.5	ab17162 - H3 P16hyd + P30hyd
B11	7.1	ab17567 - H4 K20me3
B12	4.5	ab17587 - H1.4 K25me3
B13	4.3	ab175905 - H2B K11ac
B14	4.5	ab176211 - H2B K23ac
B15	4.7	ab17632 - H4 K5bio
B16	4.6	ab1771 - H3 K9me1
B17	3.9	ab1773 - H3 K9me3
B18	6.3	ab17770 - H4 R3me1
B19	5.8	ab1781 - H3 K27me2
B20	4.6	ab1783 - H3 K36me1
B21	3.7	ab1784 - H3 K36me2
B22	4.5	ab1785 - H3 K36me3
B23	18.0	ab178667 - ScH2A unmod(1-30)
B24	23.6	ab178668 - ScH2A unmod(102-132)
B25	6.8	ab178670 - H2A.Z unmod(1-30)
C1	5.7	ab178833 - Flag
C2	3.6	ab178835 - H2B R99mutO
C3	7.7	ab178836 - H4 R23mutO
C4	5.1	ab178837 - H4 G48mutD
C5	3.8	ab178987 - H1.4 unmod(1-30)
C6	4.9	ab178988 - H1.4 S1p
C7	3.3	ab178989 - H1.4 T3p
C8	3.9	ab178990 - H1.4 S1p + T3p
C9	3.4	ab178991 - H1.4 K25me1
C10	3.0	ab178992 - H1.4 K25me3 + S26p
C11	4.0	ab178993 - H1.4 unmod(130-160)

C12	62.9	ab178994 - H2A unmod(1-30)
C13	5.3	ab178995 - H2A S1p
C14	6.5	ab178996 - H2A S1p + K5ac
C15	4.8	ab178997 - H2A S1p + R3me2a + K5ac
C16	7.4	ab178998 - H2A S1p + R3me2s + K5ac
C17	7.4	ab178999 - H2A S1p + R3cit + K5ac
C18	1542.7	ab179000 - H2A S1ac + S1p + K5ac + K9ac + K13ac + K15ac
C19	2119.5	ab179001 - H2A S1p + K5ac + K9ac + K13ac + K15ac
C20	46.8	ab179002 - H2A R3cit
C21	22.2	ab179003 - H2A R3cit + K5ac
C22	7.2	ab179004 - H2A R3me2a + K5ac
C23	6.6	ab179005 - H2A R3me2s + K5ac
C24	2423.4	ab179006 - H2A S1ac + K5ac + K9ac + K13ac + K15ac
C25	5744.5	ab179007 - H2A K5ac + K9ac + K13ac + K15ac
D1	2.8	ab179008 - H2A K9me2
D2	4.4	ab179009 - H2A K9me3
D3	3.5	ab179010 - H2A R11me1
D4	5.5	ab179011 - H2A R11me2s
D5	3.8	ab179012 - H2A R11me2a
D6	12.5	ab179013 - H2A K13ac
D7	4.0	ab179014 - H2A K15ac
D8	3.9	ab179015 - H2A R17me1
D9	3.9	ab179016 - H2A R17me2s
D10	4.0	ab179017 - H2A R17me2a
D11	10.1	ab179018 - H2A unmod(20-50)
D12	3.9	ab179019 - H2A P26hyd
D13	4.8	ab179020 - H2A R29me1
D14	4.4	ab179021 - H2A R29me2s
D15	5.5	ab179022 - H2A R29me2a
D16	5.1	ab179023 - H2A K36ac
D17	8.1	ab179024 - H2A K36cr
D18	6.5	ab179025 - H2A Y39hyd
D19	8.3	ab179026 - H2A R42me1
D20	5.3	ab179027 - H2A R42me2s
D21	6.1	ab179028 - H2A R42me2a
D22	11.4	ab179029 - H2A unmod(71-103)
D23	6.7	ab179030 - H2A R77me1
D24	380.1	ab179031 - H2A R77me2a
D25	37.1	ab179032 - H2A R88me1
E1	20.9	ab179033 - H2A R88me2s
E2	31.3	ab179034 - H2A R88me2a
E3	22.9	ab179035 - H2A K95me1
E4	21.6	ab179036 - H2A K95me3
E5	71.5	ab179037 - H2A unmod(92-123)
E6	3.6	ab179038 - H2A T101GlcNAc
E7	3.6	ab179039 - H2A Q104me1
E8	4.2	ab179040 - H2A unmod(112-129)
E9	4.3	ab179041 - H2A K118me1
E10	2.6	ab179042 - H2A K118me2
E11	4.2	ab179043 - H2A K118me3
E12	4.0	ab179044 - H2A K118fo
E13	4.1	ab179045 - H2A K118cr
E14	4.0	ab179046 - H2A K119cr
E15	5.6	ab179047 - H2A K118ubiq
E16	6.0	ab179048 - H2A K119ubiq
E17	3.3	ab179049 - H2A T120p
E18	4.1	ab179050 - H2A K125me1
E19	4.1	ab179051 - H2A K125me2
E20	8.0	ab179052 - H2A K125cr
E21	152.2	ab179053 - H2A.X unmod(1-30)
E22	296.3	ab179054 - H2A.X K5ac
E23	24.5	ab179055 - H2A.X unmod(114-142)
E24	24.5	ab179056 - H2A.X K118ubiq
E25	21.1	ab179057 - H2A.X K119ubiq

F1	16.6	ab179058 - H2A.X T120p
F2	17.5	ab179059 - H2A.X K127ubiq
F3	17.8	ab179060 - H2A.X S139p + Y142p
F4	12.3	ab179061 - H2A.Z K4ac
F5	8.8	ab179062 - H2A.Z K7ac
F6	7.9	ab179063 - H2A.Z K11ac
F7	6.5	ab179064 - H2A.Z K4ac + K7ac + K11ac
F8	7.5	ab179065 - H2B unmod(1-30)
F9	6.4	ab179066 - H2B K5ac
F10	5.9	ab179067 - H2B K5cr
F11	4.9	ab179068 - H2B K5fo
F12	5.2	ab179069 - H2B P10hyd
F13	3.2	ab179070 - H2B K11me1
F14	3.7	ab179071 - H2B K11me2
F15	3.9	ab179072 - H2B K11me3
F16	3.8	ab179073 - H2B K11cr
F17	3.9	ab179074 - H2B K12ac
F18	4.8	ab179075 - H2B K12cr
F19	3.7	ab179076 - H2B K12me1
F20	3.8	ab179077 - H2B K12me2
F21	5.4	ab179078 - H2B K12me3
F22	6.8	ab179079 - H2B K12ac + K15ac
F23	6.5	ab179080 - H2B S14p
F24	6.9	ab179081 - H2B K15ac
F25	7.7	ab179082 - H2B K15cr
G1	6.2	ab179083 - H2B K15me1
G2	6.7	ab179084 - H2B K15me2
G3	7.9	ab179085 - H2B K15me3
G4	8.3	ab179086 - H2B K16ac
G5	6.6	ab179087 - H2B K16cr
G6	7.7	ab179088 - H2B K20ac
G7	5.1	ab179089 - H2B K20cr
G8	5.9	ab179090 - H2B K20me1
G9	4.7	ab179091 - H2B K20me2
G10	5.2	ab179092 - H2B K20me3
G11	3.9	ab179093 - H2B K23cr
G12	3.8	ab179094 - H2B K23me1
G13	5.1	ab179095 - H2B K23me3
G14	13.3	ab179096 - H2B K24ac
G15	6.7	ab179097 - H2B unmod(26-56)
G16	4.0	ab179098 - H2B K34cr
G17	3.5	ab179099 - H2B K34fo
G18	3.2	ab179100 - H2B Y37hyd
G19	4.2	ab179101 - H2B Y37p
G20	3.6	ab179102 - H2B K43me2
G21	4.5	ab179103 - H2B K43me3
G22	4.7	ab179104 - H2B K46fo
G23	7.2	ab179105 - H2B K46ac
G24	7.1	ab179106 - H2B K46me1
G25	6.8	ab179107 - H2B K46me2
H1	7.0	ab179108 - H2B K46me3
H2	6.5	ab179109 - H2B unmod(51-81)
H3	7.0	ab179110 - H2B K57me1
H4	6.9	ab179111 - H2B K57me2
H5	7.3	ab179112 - H2B K57me3
H6	9.9	ab179113 - H2B unmod(73-103)
H7	4.5	ab179114 - H2B R79me1
H8	8.8	ab179115 - H2B R79me2s
H9	3.0	ab179116 - H2B R79me2a
H10	6.3	ab179117 - H2B Y83hyd
H11	8.9	ab179118 - H2B K85me1
H12	4.5	ab179119 - H2B K85me2
H13	5.4	ab179120 - H2B K85me3
H14	4.8	ab179121 - H2B unmod(94-125)

H15	5.8	ab179122 - H2B R99me1
H16	3.1	ab179123 - H2B R99me2s
H17	4.6	ab179124 - H2B R99me2a
H18	4.0	ab179125 - H2B K108fo
H19	3.7	ab179126 - H2B K116fo
H20	4.2	ab179127 - H2B K116me1
H21	2.8	ab179128 - H2B K116me2
H22	4.9	ab179129 - H2B K116me3
H23	2.8	ab179130 - H2B K116ubiq
H24	4.0	ab179131 - H2B K120ac
H25	3.2	ab179132 - H2B K120fo
I1	3.7	ab179133 - H2B K120ubiq
I2	13.5	ab179134 - H2B K125ubiq
I3	7.4	ab179135 - ScH2B unmod(102-132)
I4	10.6	ab179136 - H3 unmod(1-33)
I5	10.0	ab179137 - H3 R2cit
I6	6.7	ab179138 - H3 R2cit + K4me2
I7	7.5	ab179139 - H3 R2cit + K4me3
I8	32.8	ab179140 - H3 R2cit + K4me3 + K9ac + K14ac + K18ac
I9	6.1	ab179141 - H3 R2cit + K4ac + K9ac + K14ac + K18ac
I10	3.8	ab179142 - H3 R2me1
I11	3.8	ab179143 - H3 R2me1 + K4me2
I12	3.9	ab179144 - H3 R2me1 + K4me3
I13	6.3	ab179145 - H3 R2me1 + K4me3 + K9ac + K14ac + K18ac
I14	2.5	ab179146 - H3 R2me2s
I15	5.1	ab179147 - H3 R2me2s + K4me2
I16	7.4	ab179148 - H3 R2me2s + K4me3
I17	7.6	ab179149 - H3 R2me2s + K4me3 + K9ac + K14ac + K18ac
I18	9.6	ab179150 - H3 R2me2a
I19	3.8	ab179151 - H3 R2me2a + K4me2
I20	7.6	ab179152 - H3 R2me2a + K4me3
I21	3.6	ab179153 - H3 R2me2a + K4me3 + S10p
I22	7.6	ab179154 - H3 R2me2a + T3p + K4me3
I23	6.8	ab179155 - H3 R2me2a + K4ac + K9ac + K14ac + K18ac
I24	10.5	ab179156 - H3 R2me2a + K4me3 + K9ac + K14ac + K18ac
I25	4.7	ab179157 - H3 R2me2a + T3p + K4me3 + K9ac + K14ac + K18ac
J1	11.6	ab179158 - H3 R2me2a + K4me3 + K9ac + S10p + K14ac + K18ac
J2	4.9	ab179159 - H3 T3p
J3	6.8	ab179160 - H3 T3p + K4me1
J4	5.7	ab179161 - H3 T3p + K4me2
J5	3.4	ab179162 - H3 T3p + K4me3 + K9ac + K14ac + K18ac
J6	4.5	ab179163 - H3 K4ac
J7	5.4	ab179164 - H3 K4ac + K9ac
J8	4.8	ab179165 - H3 K4ac + K14ac
J9	9.3	ab179166 - H3 K4ac + K18ac
J10	4.5	ab179167 - H3 K4ac + K9me3
J11	5.2	ab179168 - H3 K4ac + K9ac + K14ac
J12	3.9	ab179169 - H3 K4ac + K9ac + K18ac
J13	11.4	ab179170 - H3 K4ac + K14ac + K18ac
J14	10.4	ab179171 - H3 K4ac + K9ac + K14ac + K18ac
J15	7.7	ab179172 - H3 K4ac + K9me1 + K14ac + K18ac
J16	10.9	ab179173 - H3 K4ac + K9me2 + K14ac + K18ac
J17	3.8	ab179174 - H3 K4ac + K9me3 + K14ac + K18ac
J18	8.9	ab179175 - H3 K4ac + K9ac + S10p + K14ac + K18ac
J19	5.3	ab179176 - H3 K4cr
J20	3.9	ab179177 - H3 K4me1 + K9me2
J21	5.6	ab179178 - H3 K4me1 + K9ac + K18ac
J22	7.5	ab179179 - H3 K4me1 + K9ac + K14ac + K18ac
J23	5.6	ab179180 - H3 K4me2 + T6p
J24	9.2	ab179181 - H3 K4me2 + K9me2
J25	9.6	ab179182 - H3 K4me2 + K9ac + K18ac
K1	4.4	ab179183 - H3 K4me2 + K9ac + K14ac + K18ac
K2	4.6	ab179184 - H3 K4me3 + T6p
K3	5.3	ab179185 - H3 K4me3 + K9ac

K4	3.6	ab179186 - H3 K4me3 + K9me2
K5	3.7	ab179187 - H3 K4me3 + K9me3
K6	3.7	ab179188 - H3 K4me3 + S10p
K7	4.8	ab179189 - H3 K4me3 + K14ac
K8	5.8	ab179190 - H3 K4me3 + K18ac
K9	3.8	ab179191 - H3 K4me3 + R8me2s + K9me3
K10	7.5	ab179192 - H3 K4me3 + R8me2a + K9me3
K11	5.1	ab179193 - H3 K4me3 + K9ac + S10p
K12	7.5	ab179194 - H3 K4me3 + K9ac + K14ac
K13	12.5	ab179195 - H3 K4me3 + K9ac + K18ac
K14	19.1	ab179196 - H3 K4me3 + K14ac + K18ac
K15	6.9	ab179197 - H3 K4me3 + K9ac + K14ac + K18ac
K16	7.8	ab179198 - H3 K4me3 + K9ac + S10p + K14ac + K18ac
K17	6.5	ab179199 - H3 K4me3 + T6p + K9ac + K14ac + K18ac
K18	8.4	ab179200 - H3 K4me2 + T6p + K9ac + K14ac + K18ac
K19	77.7	ab179201 - H3 K4bio + K9bio + K18bio
K20	4.6	ab179202 - H3 T6p + K9me3
K21	5.4	ab179203 - H3 T6p + R8me2a + K9me3
K22	8.7	ab179204 - H3 T6p + K9ac + K14ac + K18ac
K23	6.9	ab179205 - H3 R8me1
K24	5.7	ab179206 - H3 R8me1 + K9me1
K25	5.3	ab179207 - H3 R8me1 + K9me2
L1	3.4	ab179208 - H3 R8me1 + K9me3
L2	5.5	ab179209 - H3 R8me2s
L3	3.5	ab179210 - H3 R8me2s + K9me1
L4	4.1	ab179211 - H3 R8me2s + K9me2
L5	3.8	ab179212 - H3 R8me2s + K9me3
L6	6.6	ab179213 - H3 R8me2a
L7	4.1	ab179214 - H3 R8me2a + K9me1
L8	7.2	ab179215 - H3 R8me2a + K9me2
L9	3.9	ab179216 - H3 R8me2a + K9me3
L10	7.6	ab179217 - H3 R8cit
L11	6.5	ab179218 - H3 K9ac + K14ac
L12	15.8	ab179219 - H3 K9ac + K18ac
L13	6.7	ab179220 - H3 K9ac + K14ac + K18ac
L14	6.4	ab179221 - H3 K9ac + K14ac + K18ac + K23ac + K27ac
L15	3.7	ab179222 - H3 K9ac + S10p
L16	4.0	ab179223 - H3 K9cr
L17	4.1	ab179224 - H3 K9me2
L18	5.1	ab179225 - H3 K9me2 + S10p
L19	5.8	ab179226 - H3 K9me2 + K27me2
L20	4.0	ab179227 - H3 S10p + S28p
L21	3.4	ab179228 - H3 S10p K14ac
L22	4.8	ab179229 - H3 K14me1
L23	3.9	ab179230 - H3 K14me2
L24	3.0	ab179231 - H3 K14me3
L25	3.5	ab179232 - H3 K14ac + K18ac
M1	5.0	ab179233 - H3 R17cit
M2	3.6	ab179234 - H3 R17me1
M3	5.6	ab179235 - H3 R17me2s
M4	6.2	ab179236 - H3 R17me2a + K18ac
M5	5.4	ab179237 - H3 K18cr
M6	5.3	ab179238 - H3 K18me1
M7	6.7	ab179239 - H3 K18me2
M8	10.8	ab179240 - H3 K18me3 + K36me3
M9	3.5	ab179241 - H3 K23cr
M10	5.6	ab179242 - H3 K23me1
M11	10.0	ab179243 - H3 K23me2
M12	7.8	ab179244 - H3 K23me3
M13	5.3	ab179245 - H3 T22p + K23me3
M14	4.7	ab179246 - H3 K23pr
M15	5.1	ab179247 - H3 unmod(20-51)
M16	4.3	ab179248 - H3 R26me1
M17	3.5	ab179249 - H3 R26me2a

M18	3.2	ab179250 - H3 R26me2a + K27me1
M19	2.7	ab179251 - H3 R26me2a + K27me3
M20	6.4	ab179252 - H3 K27ac + S28p
M21	5.7	ab179253 - H3 K27cr
M22	5.3	ab179254 - H3 K27me1
M23	5.2	ab179255 - H3 K27me2 + S28p
M24	6.0	ab179256 - H3 K27me3
M25	5.3	ab179257 - H3 K27me3 + S28p
N1	4.5	ab179258 - H3 P30hpro
N2	5.6	ab179259 - H3 P30(5,5me2)
N3	3.5	ab179260 - H3 P30
N4	3.9	ab179261 - H3 P30mutA
N5	3.5	ab179262 - H3 T32GlcNAc
N6	7.9	ab179263 - H3 K36ac
N7	7.3	ab179264 - H3 P38hpro
N8	4.5	ab179265 - H3 P38(5,5me2)
N9	5.6	ab179266 - H3 P38
N10	4.8	ab179267 - H3 P38mutA
N11	7.5	ab179268 - H3 Y41p
N12	7.9	ab179269 - H3 unmod(50-80)
N13	5.1	ab179270 - H3 K56ac
N14	5.4	ab179271 - H3 K56cr
N15	4.3	ab179272 - H3 K56fo
N16	5.4	ab179273 - H3 K56me1
N17	9.5	ab179274 - H3 K56me2
N18	4.9	ab179275 - H3 K56me3
N19	3.9	ab179276 - H3 R63me1
N20	4.3	ab179277 - H3 R63me2s
N21	6.2	ab179278 - H3 R63me2a
N22	3.0	ab179279 - H3 K64ac
N23	5.6	ab179280 - H3 K64me2
N24	6.4	ab179281 - H3 K64me3
N25	4.4	ab179282 - H3 unmod(73-103)
O1	683.1	ab179283 - H3 K79ac
O2	26.7	ab179284 - H3 K79fo
O3	5.5	ab179285 - H3 T80p
O4	5.7	ab179286 - H3 K79me3 + T80p
O5	4.9	ab179287 - H3 unmod(99-135)
O6	6.4	ab179288 - H3 T118O-ac
O7	5.4	ab179289 - H3 K122fo
O8	5.8	ab179290 - H3 K122me1
O9	6.7	ab179291 - H3 R128me1
O10	6.0	ab179292 - H3 R128me2s
O11	12.2	ab179293 - H3.3 unmod(10-40)
O12	4.8	ab179294 - H3.3 S28p
O13	5.8	ab179295 - H3.3 S28p + S31p
O14	4.8	ab179296 - H3.3 S31p
O15	9.4	ab179297 - H4 unmod(1-30)
O16	4.4	ab179298 - H4 S1p + R3me1
O17	5.1	ab179299 - H4 S1p + R3me2s
O18	6.6	ab179300 - H4 S1p + R3me2a
O19	17.5	ab179301 - H4 S1p + K5ac + K8ac + K12ac + K16ac
O20	11.9	ab179302 - H4 R3cit
O21	14.0	ab179303 - H4 R3me1 + K5ac + K8ac + K12ac + K16ac + K20ac
O22	14.2	ab179304 - H4 R3me2s
O23	12.7	ab179305 - H4 R3me2s + K5ac + K8ac + K12ac + K16ac + K20ac
O24	32.1	ab179306 - H4 R3me2a
O25	5.4	ab179307 - H4 R3me2a + K5ac
P1	89.6	ab179308 - H4 R3me2a + K5ac + K8ac + K12ac + K16ac + K20ac
P2	6.5	ab179309 - H4 K5ac
P3	7.3	ab179310 - H4 K5ac + K8ac
P4	6.6	ab179311 - H4 K5ac + K12ac
P5	16.5	ab179312 - H4 K5ac + K16ac
P6	8.0	ab179313 - H4 K5ac + K8ac + K12ac

P7	9.3	ab179314 - H4 K5ac + K8me1 + K12ac
P8	14.5	ab179315 - H4 K5ac + K8ac + K12ac + K16ac
P9	9.2	ab179316 - H4 K5ac + K8me1 + K12ac + K16ac
P10	35.3	ab179317 - H4 K5ac + K8ac + K12me1 + K16ac
P11	20.3	ab179318 - H4 K5me1
P12	5.3	ab179319 - H4 K5me1 + K8ac + K12me1
P13	7.7	ab179320 - H4 K5me1 + K8me1 + K12me1
P14	5.5	ab179321 - H4 K5me1 + K8ac + K12ac + K16ac
P15	26.1	ab179322 - H4 K5me2
P16	8.0	ab179323 - H4 K5me3
P17	45.4	ab179324 - H4 K5cr
P18	6.3	ab179325 - H4 K8ac + K12ac
P19	9.3	ab179326 - H4 K8ac + K16ac
P20	7.2	ab179327 - H4 K8cr
P21	7.3	ab179328 - H4 K8me1
P22	5.4	ab179329 - H4 K12ac
P23	5.7	ab179330 - H4 K12ac + K16ac
P24	11.0	ab179331 - H4 K12ac + K16ac + K20me2
P25	6.2	ab179332 - H4 K12ac + K16ac + K20me3
Q1	199.2	ab179333 - H4 K5mutQ + K8mutQ + K12ac + K16mutQ + K20mutQ
Q2	7.2	ab179334 - H4 K12cr
Q3	6.8	ab179335 - H4 K12fo
Q4	5.3	ab179336 - H4 K12me1
Q5	8.8	ab179337 - H4 K12me2
Q6	5.4	ab179338 - H4 K12me3
Q7	8.3	ab179339 - H4 K16ac
Q8	5.9	ab179340 - H4 K16ac + K20ac
Q9	6.3	ab179341 - H4 K16cr
Q10	5.9	ab179342 - H4 K16me1
Q11	5.7	ab179343 - H4 K16me3
Q12	10.8	ab179344 - H4 R17me1
Q13	9.8	ab179345 - H4 R17me2s
Q14	23.3	ab179346 - H4 R17me2a
Q15	10.4	ab179347 - H4 R19me1
Q16	6.3	ab179348 - H4 R19me2s
Q17	7.8	ab179349 - H4 R19me2a
Q18	6.5	ab179350 - H4 K20ac
Q19	9.5	ab179351 - H4 R23me1
Q20	6.7	ab179352 - H4 R23me2s
Q21	6.7	ab179353 - H4 R23me2a
Q22	5.3	ab179354 - H4 unmod(24-54)
Q23	5.1	ab179355 - H4 K31fo
Q24	5.2	ab179356 - H4 K31me1
Q25	8.0	ab179357 - H4 K31me2
R1	5.0	ab179358 - H4 K31me3
R2	6.1	ab179359 - H4 P32hyd
R3	5.4	ab179360 - H4 R35me1
R4	7.5	ab179361 - H4 R35me2s
R5	6.3	ab179362 - H4 R35me2a
R6	7.6	ab179363 - H4 S47O-ac
R7	5.6	ab179364 - H4 S47GlcNAc
R8	7.1	ab179365 - H4 S47p
R9	7.0	ab179366 - H4 unmod(45-75)
R10	8.6	ab179367 - H4 Y51hyd
R11	7.3	ab179368 - H4 R55me1
R12	4.4	ab179369 - H4 R55me2s
R13	5.3	ab179370 - H4 R55me2a
R14	4.9	ab179371 - H4 K59fo
R15	5.2	ab179372 - H4 K59me1
R16	5.6	ab179373 - H4 K59me2
R17	6.2	ab179374 - H4 R67me1
R18	4.4	ab179375 - H4 R67me2s
R19	5.9	ab179376 - H4 R67me2a
R20	4.3	ab179377 - H4 unmod(71-102)

R21	4.8	ab179378 - H4 K77me1
R22	5.3	ab179379 - H4 K77me2
R23	5.8	ab179380 - H4 K77me3
R24	7.0	ab179381 - H4 K79fo
R25	6.7	ab179382 - H4 K79me1
S1	4.7	ab179383 - H4 K79me3
S2	5.8	ab179384 - H4 T82O-ac
S3	7.1	ab179385 - H4 Y88hyd
S4	6.4	ab179386 - H4 Y88nitro
S5	5.2	ab179387 - H4 K91fo
S6	25.4	ab179403 - H2A K5ac
S7	812.0	ab179404 - H2A K9me1
S8	4.4	ab18504 - H2B S32p
S9	3.4	ab19828 - ScH2A S129p
S10	4.5	ab20631 - H3 R26cit
S11	6.5	ab20635 - H3 T3p + K4me3
S12	2.8	ab21043 - H3 K37me1
S13	2.7	ab21998 - H1.4 K25me2
S14	5.0	ab22079 - H4 K59me3
S15	5.7	ab22080 - H3 K122me2
S16	7.7	ab22082 - H3 K122me3
S17	5.2	ab22214 - H3 K37me2
S18	17.7	ab22398 - H2A R3me2a
S19	206.5	ab22399 - H2A R3me2s
S20	10.8	ab23539 - H3 T22p
S21	4.5	ab24003 - H3 K18ac
S22	7.7	ab24404 - H3 K27ac
S23	3.9	ab24417 - H3 K37me3
S24	3.6	ab24444 - H3 T11p
S25	2.7	ab24659 - H2A K127me3
T1	3.5	ab24660 - H2A K125me3
T2	5.1	ab26307 - H3 T45p
T3	3.7	ab26342 - H2B K5me2
T4	6.6	ab26739 - H4 T30p
T5	5.9	ab27896 - H4 T96p
T6	8.0	ab27898 - H4 T80p
T7	224.9	ab28756 - H2A R3me1
T8	9.7	ab32876 - H3 R2cit + R8cit + R17cit
T9	8.2	ab33047 - H2A K127me2
T10	7.5	ab33505 - H3 T118p
T11	6.3	ab33594 - H3 R128me2a
T12	5.6	ab34466 - H3 K122ac
T13	3.5	ab36927 - H2A S122p
T14	4.3	ab4555 - H3 K79me1
T15	4.9	ab4556 - H3 K79me2
T16	4.2	ab4557 - H3 K79me3
T17	5.9	ab4560 - H4 K79me2
T18	4.4	ab46854 - H3 K18me3
T19	21.5	ab48359 - H3 K23ac
T20	6.7	ab53244 - H3 K64me1
T21	2857.2	ab54016 - H2A K9ac
T22	4.0	ab5499 - H3 S28p
T23	6.6	ab7768 - H3 K4me2
T24	3.3	ab89018 - H1.4 S26p
T25	5.8	ab92374 - H3 K4me3 (HEX1)
U1	6.6	ab95444 - H4 K16me2