





Anti-Histone H2B (acetyl K16) antibody [EPR17598] – ChIP Grade (ab177427)

Target	Histone H2B (acetyl K16)
Sample	HeLa
Treatment	None
Antibody amount	4 µg
Cell number	1×10 ⁷
Sequencing	NovaSeq 6000 (PE 150)
Peak calling	MACS 2.1.0
Reference genome	hg38








Chromosome 16 (active region)





Copyright (c) 2021 Abcam plc.



Copyright (c) 2021 Abcam plc.



Copyright (c) 2021 Abcam plc.



Copyright (c) 2021 Abcam plc.