


Anti-ARID1A antibody [EPR13501-73] - (ab182561)


Application	ChIC/CUT&RUN-seq
Target	ARID1A
Sample	HCT116
Treatment	None
Antibody amount	5µg
Cell number	2.5X10 ⁵
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	hg38
Reads	10 million

Chromosome 4




Copyright (c) 2023 Abcam plc.

Chromosome 20




Copyright (c) 2023 Abcam plc.

Chromosome 17



Copyright (c) 2023 Abcam plc.

Chromosome 3



Copyright (c) 2023 Abcam plc.