


## Anti-CTCF antibody [EPR18253] – ChIP Grade (ab188408)


|                  |                      |
|------------------|----------------------|
| Application      | ChIC/CUT&RUN-seq     |
| Target           | CTCF                 |
| Sample           | F9                   |
| Treatment        | None                 |
| Antibody amount  | 5µg                  |
| Cell number      | 2.5X10 <sup>5</sup>  |
| Sequencing       | NovaSeq 6000 (PE150) |
| Peak calling     | MACS 2.1.0           |
| Reference genome | mm39                 |
| Reads            | 10 million           |

### Chromosome 10




Copyright (c) 2023 Abcam plc.

### Chromosome 12




Copyright (c) 2023 Abcam plc.

### Chromosome X



Copyright (c) 2023 Abcam plc.

### Chromosome 6



Copyright (c) 2023 Abcam plc.