


Anti-ATF3 antibody [EPR19488] - (ab207434)


Application	ChIC/CUT&RUN-seq
Target	ATF3
Sample	Human ATF3 knockout HeLa cell line (ab264908) Human wild-type HeLa cell line (ab255448)
Treatment	None
Antibody amount	5µg
Cell number	2.5X10 ⁵
Antibody used for Assay Quality Control	Anti-CTCF (ab188408)
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	hg38
Reads	10 million

Chromosome 12 zoom out




Copyright (c) 2023 Abcam plc.

ab207434 Chromosome 1




Copyright (c) 2023 Abcam plc.

ab207434 Chromosome 6



Copyright (c) 2023 Abcam plc.

ab207434 Chromosome 1



ab207434 Chromosome 3

