


Anti- Smad2+Smad3 antibody [EPR19557] - (ab207447)


Application	ChIC/CUT&RUN-seq
Target	Smad2+Smad3
Sample	HaCaT
Treatment	TGF- β 1(7ng/ml 1h)
Antibody amount	5 μ g
Cell number	2.5X10 ⁵
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	hg38
Reads	10 million

Chromosome 19




Copyright (c) 2023 Abcam plc.

Chromosome 20




Copyright (c) 2023 Abcam plc.

Chromosome 15



Copyright (c) 2023 Abcam plc.

Chromosome 5



Copyright (c) 2023 Abcam plc.