



Anti-Nanog antibody [EPR20694] – ChIP Grade (ab214549)

Application	ChIC/CUT&RUN-seq
Target	Nanog
Sample	F9
Treatment	None
Antibody amount	5µg
Cell number	2.5X10 ⁵
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	Mm39
Reads	10 million


Chromosome 12



Chromosome 6




Chromosome 3



Copyright (c) 2023 Abcam plc.

Chromosome X



Copyright (c) 2023 Abcam plc.