



Anti-SP1 antibody [EPR22648-50] - ChIP Grade (ab231778)

Application	ChIC/CUT&RUN-seq
Target	SP1
Knock Out (KO) cell line	Human SP1 knockout HeLa cell line (ab265519)
Wild Type (WT) cell line	wild-type HeLa cell line (ab255928)
Treatment	None
Antibody amount	5µg
Cell number	2.5X10 ⁵
Antibody used for Assay Quality Control	Anti-CTCF(ab188408)
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	hg38
Reads	10 million

Chromosome 12 zoom out




ab231778 Chromosome 12



Copyright (c) 2022 Abcam plc.

ab231778 Chromosome 14



Copyright (c) 2022 Abcam plc.

ab231778 Chromosome 7

