




Anti-Histone H3 (acetyl K64) antibody [EPR20713] - ChIP Grade (ab214808)

Target	Histone H3 (acetyl K64)
Sample	HeLa
Treatment	None
Antibody amount	4 µg
Cell number	1×10 ⁷
Sequencing	NovaSeq 6000 (PE 150)
Peak calling	MACS 2.1.0
Reference genome	hg38



30,000 kb 30,025 kb 30,050 kb 30,075 kb 30,100 kb 30,125 kb

Chromosome 16 (active region)



Copyright (c) 2021 Abcam plc.


17,690 kb 17,700 kb 17,710 kb 17,720 kb 17,730 kb 17,740 kb 17,750 kb

Chromosome 11 (inactive region)



Copyright (c) 2021 Abcam plc.

73,300 kb 73,400 kb 73,500 kb 73,600 kb



Copyright (c) 2021 Abcam plc.

Copyright (c) 2021 Abcam plc.