



Anti- Smad2 (phospho T8) + Smad3 (phospho T8) antibody [EPR23682-64] - (ab254407)

Application	ChIC/CUT&RUN-seq
Target	Smad2 (phospho T8) + Smad3 (phospho T8)
Sample	HaCaT
Treatment	None
Antibody amount	5µg
Cell number	2.5X10^5
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	hg38
Reads	10 million


Chromosome 20



Chromosome 5




Chromosome 6



Copyright (c) 2023 Abcam plc.

Chromosome 15



Copyright (c) 2023 Abcam plc.