


Anti- KAT3B/P300 antibody [EPR23753-55] - (ab259330)


Application	ChIC/CUT&RUN-seq
Target	KAT3B/P300
Sample	HeLa
Treatment	None
Antibody amount	5µg
Cell number	2.5X10 ⁵
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	hg38
Reads	10 million

Chromosome 5




Copyright (c) 2023 Abcam plc.

Chromosome 1




Copyright (c) 2023 Abcam plc.

Chromosome 17




Copyright (c) 2023 Abcam plc.

Chromosome 5



Copyright (c) 2023 Abcam plc.

Chromosome 6



Copyright (c) 2023 Abcam plc.