


Anti- RUNX1 / AML1 antibody [EPR23309-113] - ChIP Grade (ab272456)


Application	ChIC/CUT&RUN-seq
Target	RUNX1 / AML1
Sample	K-562
Treatment	None
Antibody amount	5µg
Cell number	2.5X10 ⁵
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	hg38
Reads	10 million

Chromosome 3




Copyright (c) 2023 Abcam plc.

Chromosome 10




Copyright (c) 2023 Abcam plc.

Chromosome 11



Copyright (c) 2023 Abcam plc.

Chromosome 1



Copyright (c) 2023 Abcam plc.